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MEGBI Training Course III

Vaccine production techniques: Egg based and Cell based amplification of influenza virus

1 day training course for researchers and technical staff

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Authors:

Noha Abdulwahab, Ghina al-Eter, Omran Zaki, Layal Chbib, and Mirna Khoder

Institute for Genetic Engineering, Ecology and Health (IGEEH)

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**Middle East Genetics and Biotechnology
Institute (MEGBI)**

Ras Nhache, Batroun, Lebanon,

www.aecenar.com/institutes/megbi

Email: info@aecenar.com

Remark: All photos in part II (Cell based virus propagation) are taken from laboratory work in MEGBI, Ras Nhache, Lebanon.

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I. Egg Based Influenza Virus Propagation

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Chickens are susceptible to many infectious diseases. One of the most important of these is the viral disease known as influenza which is caused for examples by the actual H1N1 virus. For this reasons viruses can be propagated will in chicken eggs for vaccine production purpose.

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Eggs for our work can be purchase from normal chicken farms.the eggs must be 9-10 days old when purchased.

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Then the virus probe is inoculated into the eggs where the virus is propagated while the eggs are incubated. The work must be done under very clean circumstances and atmosphere to prevent contamination.

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Testing the actual state of the eggs is done by candling.

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Laboratory staff should be familiar with and have practiced the following skills prior to the commencement of H1N1 disease vaccine production - this manual does not contain further details about these skills:

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- Aseptic technique.
- Sterilization by autoclaving and hot air of glassware and discarded materials.

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An order can be placed for the delivery of the eggs. It is useful if the person responsible for placing the orders and receiving the eggs keeps records. The following information should be recorded in a notebook set aside for this purpose.

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Appropriate chemical disinfectants must be used for cleaning equipment, materials and work surfaces. All laboratory wastes must be assigned to a category and placed in clearly labeled

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bins from where the waste will be disposed of appropriately.

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Alcohol

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A 70% volume/volume (v/v) solution of alcohol diluted with water is useful for wiping down benches and disinfecting the outside of eggs before inoculation and harvesting of allantoic fluid. The addition of 2 percent iodine will increase the effectiveness of this solution.

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Note that 70 percent solution of alcohol is flammable!

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Chlorine

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There are several chlorine compounds that are used as disinfectants. Sodium hypochlorite (NaOCl) or household bleach is readily available and cheap.

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Soaking overnight in a 2 percent solution of chlorine is useful for disinfecting plastic materials. Note that commercial bleach contains 12 to 14 percent hypochlorite when manufactured but this concentration deteriorates with time.

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Note that chlorine damages fabric and corrodes many metals!

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Always read the instructions before using disinfectants and cleaning reagents!

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Many vaccine production centres will already have large commercial incubators installed. Smaller incubators are available and are suitable for the small-scale production of vaccine.

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Inoculated eggs contain virus and should be placed in a different incubator.

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1.6 Cleaning and decontamination of incubators / تنظيف و تطهير الحضانة

Keep surfaces clean by wiping out with a wet cloth and disinfecting with 70 percent alcohol solution or a non-corrosive disinfectant.

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
Candling is the process of holding a strong light above or below the egg to observe the embryo. A candling lamp consists of a strong electric bulb covered by a plastic or aluminium container that has a handle and an aperture. The egg is placed against this aperture and illuminated by the light. If you do not have a candling lamp, improvise. Try using a torch. Candling is done in a darkened room or in an area shielded by curtains.

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Note: In some eggs the air sac will have not developed on the blunt end but half way down the egg. These eggs are not suitable for vaccine production.

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The most convenient method of propagating of influenza virus in the laboratory is by the inoculation of the allantoic cavity of embryonated eggs. All strains of virus will grow in the cell of the allantoic cavity. The virus enters these cells where it multiplies. As the cells are disrupted the virus is shed into the allantoic fluid. Virulent strains of the virus will invade cells beyond the lining of the allantoic cavity and kill the embryo. The time taken for this to occur is the basis of the "Mean Death Time Assays", which indicates the level of virulence. The avirulent strain of influenza virus will not kill embryos inoculated into the allantoic cavity.

Inoculation of the allantoic cavity of embryonated eggs is a technique used in the following procedures:

1. influenza virus vaccine production
2. Establishing the infectivity titre of a suspension of virus.
3. Isolation of virus from field specimens for laboratory Diagnosis

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¹ Taken from Food and Agriculture Organization of the United Nations (www.fao.org) document <ftp://ftp.fao.org/docrep/fao/005/ac802e/ac802e01.pdf>.

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Materials

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Eggs 9-day old or 10-day old embryonated eggs. Candle the eggs and mark the inoculation sites . Eggs should be placed in an egg rack with the inoculation site uppermost.

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- Egg shell punch or forceps.
- Cotton.
- A 70 % alcohol solution in water.
- Syringe 1 mL.
- Needles preferably 25 gauge, 16 mm.
- sticky tape or melted wax to seal the inoculation site.
- Inoculum. This must be free of microbial contamination.
- Discard tray.

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Method

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1. Use cotton wool and 70 percent alcohol to swab the end of the eggs to be inoculated. Allow the alcohol to evaporate.

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2. Swab the eggshell punch with 70 percent alcohol solution. Place used cotton wool in discard tray.

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3. Pierce a hole in the end of the egg at the marked inoculation site.

4. Attach needle to 1 mL syringe.

5. Draw inoculum into 1 mL syringe.

6. Keeping the needle and syringe vertical, place the needle through the hole in the eggshell. The needle will need to penetrate approximately 16mm into the egg to reach the allantoic cavity.

7. Inject 0.1 mL of inoculum into the egg.

8. Withdraw the needle from the egg.

9. Seal the hole in the shell with stationery tape or melted wax.

10. Discard the used needles and syringes.

11. Place the inoculated eggs into a second incubator. Check the temperature and humidity of incubator.


Figure 2: Inoculation of the allantoic cavity

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The following method describes harvesting a small sample of allantoic fluid for testing for the presence of haemagglutinin using the rapid or micro tests.

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Materials

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- Forceps or a small pair of scissors
- Absolute alcohol for flaming forceps
- Cotton.
- 70 percent alcohol solution
- Discard tray
- 50 µL micropipette and tips, a wire loop or sterile Pasteur pipettes

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1. Chill eggs at 4°C for at least two hours to kill the embryo and to reduce the contamination of the allantoic fluid with blood during harvesting.

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2. Remove sticky tape (if used to seal the eggs) and swab each egg with cotton wool soaked with 70 percent alcohol to disinfect and remove

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3. Dip the forceps or scissors in absolute alcohol and flame to sterilize. Remove the eggshell above the air space.

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4. Discard embryos that are visibly contaminated.

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5. Remove a sample of allantoic fluid from each egg. Use a micropipette and sterile tip, sterile glass pipette or a flamed loop and dispense the sample according to method being used for the test

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This is the result of the haemagglutinin part of the haemagglutinin/neuraminidase viral protein binding to receptors on the membrane of red blood cells. The linking together of the red blood cells by the viral particles results in clumping. This clumping is known as haemagglutination. Haem-agglutination is visible macroscopically and is the basis of haemagglutination tests to detect the presence of viral particles. The test does not discriminate between viral particles that are infectious and particles that are degraded and no longer able to infect cells. Both can cause the agglutination of red blood cells.

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Every time a haemagglutination test is carried out, it is necessary to test the settling pattern of the suspension of red blood cells. This involves mixing diluent with red blood cells and allowing the cells to settle.

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2. Add red blood cells and mix by gently shaking.

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3. Allow the red blood cells to settle and observe the pattern.

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4. Observe if the cells have a normal settling pattern and there is no auto-agglutination. This will be a distinct button of cells in the micro test and an even suspension with no signs of clumping in the rapid test.

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Note: The diluent used for haemagglutination tests in this manual is PBS. There should be no signs of haemolysis in the red blood cell suspension. If there are signs of haemolysis, a fresh suspension must be prepared. There should not be any sign of auto-agglutination in the red blood cell control. If an agglutination pattern is observed, discard the suspension of red blood cells. Prepare a fresh suspension and test again.

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Tissue culture is often a generic term that refers to both organ culture and cell culture and the terms are often used interchangeably. Cell cultures are derived from either primary tissue explants or cell suspensions. Primary cell cultures typically will have a finite life span in culture whereas continuous cell lines are, by definition, abnormal and are often transformed cell lines.

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When cells are taken freshly from animal tissue and placed in culture, the cultures consist of a wide variety of cell types, most of which are capable of very limited growth in vitro, usually fewer than ten divisions. These cells retain their diploid karyotype, i.e., they have the chromosome number and morphology of their tissues of origin. They also retain some of the differentiated characteristics that they possessed in vivo. Because of this, these cells support the replication of a wide range of viruses. Primary cultures derived from monkey kidneys, mouse fetuses, and chick embryos are commonly used for diagnostic purposes and

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Some primary cells can be passed through secondary and several subsequent subcultures while retaining their original morphological characteristics and karyotype. Subcultures will have fewer cell types than primary cultures. After 20 to 50 passages in vitro, these diploid cell strains usually undergo a crisis in which their growth rate slows and they eventually die out. Diploid strains of fibroblasts derived from human fetal tissue are widely used in diagnostic virology and vaccine production.

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Certain cultured cells, notably mouse fetal fibroblasts, kidney cells from various mammalian species, and human carcinoma cells, are able to survive the growth crisis and undergo indefinite propagation in vitro. After several passages, the growth rate of the culture slows down; then isolated colonies of cells begin to grow more rapidly than diploid cells, their karyotype becomes abnormal , their morphology changes, and other poorly understood changes take place that make the cells immortal. The cells are now "dedifferentiated," having lost the specialized morphology and

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and allow movement for growth. The most convenient vessels are specially-treated polystyrene plastic that are supplied sterile and are disposable. These include petri dishes, multi-well plates, microtiter plates, roller bottles, and screwcap flasks .

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3.3 Preservation and storage of tissue cells / الحفظ و التخزين

(If available) liquid N₂ is used to preserve tissue culture cells, either in the liquid phase (-196°C) or in the vapor phase (-156°C). Freezing can be lethal to cells due to the effects of damage by ice crystals, alterations in the concentration of electrolytes, dehydration, and changes in pH. To minimize the effects of freezing, several precautions are taken. First, a cryoprotective agent which lowers the freezing point, such as glycerol or DMSO, is added. A typical freezing medium is 90% serum, 10% DMSO. In addition, it is best to use healthy cells that are growing in log phase and to replace the medium 24 hours before freezing. Also, the cells are slowly cooled from room temperature to -80°C to allow the water to move out of the cells before it freezes. The optimal rate of cooling is 1°-3°C per minute.

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The Mr. Frosty is filled with 200 ml of isopropanol at room temperature and the freezing vials containing the cells are placed in the container and the container is placed in the -80°C freezer. The effect of the isopropanol is to

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allow the tubes to come to the temperature of the freezer slowly, at about 1°C per minute.

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MAINTENANCE

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Cultures should be examined daily, observing the morphology, the color of the medium and the density of the cells. A tissue culture log should be maintained that is separate from your regular laboratory notebook. The log should contain: the name of the cell line, the medium components and any alterations to the standard medium, the dates on which the cells were split and/or fed, a calculation of the doubling time of the culture (this should be done at least once during the semester), etc.

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3.4 Harvesting and refeeding culture cells / الحصاد

Cells are harvested when the cells have reached a high density which suppresses growth.

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Suspension cultures are fed by dilution into fresh medium.

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Adherent cultures that do not need to be divided can simply be fed by removing the old medium and replacing it with fresh

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When the cells become semi-confluent, several methods are used to remove the cells from the growing surface so that they can be diluted.

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A rubber spatula can be used to physically remove the cells from the growth surface. This method is quick and easy but is also disruptive to the cells and may result in significant cell death. This method is best when harvesting many different samples of cells for preparing extracts.

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Trypsin, collagenase, or pronase, usually in combination with EDTA, causes cells to detach from the growth surface. This method is fast and reliable but can damage the cell surface by digesting exposed cell surface proteins. The proteolysis reaction can be quickly terminated by the addition of complete medium containing serum

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EDTA alone can also be used to detach cells and seems to be gentler on the cells than trypsin.

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The standard procedure for detaching adherent cells is as follows:

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The viability of cells can be observed visually using an inverted phase contrast microscope. Live cells are phase bright; suspension cells are typically rounded and somewhat symmetrical; adherent cells will form projections when they attach to the growth surface. Viability can also be assessed using the vital dye, trypan blue, which is excluded by live cells but accumulates in dead cells. Cell numbers are determined using a Hemacytometer.

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Assume all cultures are hazardous since they may harbor latent viruses or other organisms that are uncharacterized.

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The following safety precautions should also be observed:

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The cells should be monitored daily for morphology and growth characteristics, fed every 2 to 3 days, and subcultured when necessary. A minimum of two 25 cm² flasks should be carried for each cell line. Each time the cells are subcultured, a viable cell count should be done, the subculture dilutions should be noted, and, after several passages, a doubling time determined. As soon as you have enough cells, several vials should be frozen away and stored in liquid N₂. One vial from each freeze down should be thawed 1-2 weeks after freezing to check for viability. These frozen stocks will prove to be vital if any of your cultures become contaminated.

the minimal nutritional requirements of cultured cells :mixture of salts, amino acids, vitamins and cofactors, carbohydrates, and horse serum. By eliminating one component at a time , then determined which nutrients were essential for cell growth. His minimum essential medium (MEM) contains 13 amino acids (human tissue in vivo requires only eight), 8 vitamins and cofactors, glucose as an energy source, and a physiological

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salt solution that is isotonic to the cell. The pH is maintained at 7.2 to 7.4 by NaHCO_3 in equilibrium with CO_2 . The pH indicator phenol red is usually incorporated into the medium; it turns red-purple if the medium is basic, yellow if the medium is acidic, and remains red-orange if the pH is in the right range.

Serum in concentrations of 1 to 10% must be added to the medium to provide the cells with additional poorly defined factors, without which most cells will not grow. Most mammalian cells are incubated at 37°C ; avian, reptilian, and arthropod cells may grow best at higher or lower temperatures.

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When adherent cells become semi-confluent, subculture using 2 mM EDTA or trypsin/EDTA.

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VIABLE CELL NUMBERS.

Blue is one of several stains recommended for use in dye exclusion procedures for viable cell counting. This method is based on the principle that live cells do not take up certain dyes, whereas dead cells do.

- Prepare a cell suspension, either directly from a cell culture or from a concentrated or diluted suspension and combine 20 µl of cells with 20 µl of trypan blue suspension (0.4%). Mix thoroughly and allow to stand for 5-15 minutes.
- With the cover slip in place, transfer a small amount of trypan blue-cell suspension to both chambers of the Hemacytometer by carefully touching the edge of the cover slip with the pipette tip and allowing each chamber to fill by capillary action. Do not overfill or underfill the chambers.
- Starting with 1 chamber of the Hemacytometer, count all the cells in the 1 mm center square and four 1 mm corner square. Keep a separate count of viable and non-viable cells.
- If there are too many or too few cells to count, repeat the procedure either concentrating or diluting the original suspension as appropriate.
- Include cells on top and left touching middle line. Do not count cells touching middle line at bottom and right. Count 4 corner squares and middle square in both chambers

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- Each large square of the Hemacytometer, with cover-slip in place, represents a total volume of 0.1 mm^3 or 10^{-4} cm^3 . Since 1 cm^3 is equivalent to approximately 1 ml, the total number of cells per ml will be determined using the following calculations: Cells/ml = average cell count per square x dilution factor x 10^4 ;

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- Total cells = cells/ml x the original volume of fluid from which the cell sample was removed; % Cell viability = total viable cells /total cells x 100.

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Add distilled H₂O to 1 litre. Filter sterilize. Saline A is usually prepared as a 10X solution and stored at -20 °C

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
Cell based Influenza Virus propagation for vaccine production

L-cystine.HCl	0.0313
L-histidine.HCl.H ₂ O	0.042
L-isoleucine	0.052
L-leucine	0.052
L-lysine.HCl	0.0725
L-methionine	0.015
L-phenylalanine	0.032
L-threonine	0.048
L-tryptophan	0.010
L-tyrosine 2Na.2H ₂ O	0.0519
L-valine	0.001
Folic acid	0.001
myo-Inositol	0.002
Niacinamide	0.001
D-Pantothenic acid (calcium)	0.001
Pyridoxal.HCl	0.001
Riboflavin	0.0001
Thiamine.HCl	0.001
Glucose	1.0
Phenol red	0.011
NaHCO ₃	2.2

Add distilled H₂O to 1 liter. Filter sterilize. Minimum essential medium Eagle is usually purchased as a preweighed mixture or as a sterile solution. The above formulation has Earle's salts.

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- Incubator
- Centrifuge
- Forceps and scissors
- Alcohol

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
4.3 Protocol

1-Disinfect the surface of the egg over the air sac with scissors or forceps, break the shell.

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2-Sterilize forceps by dipping in alcohol and flaming. Cool forceps, then peel away the shell over the air sac

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3-Sterilize forceps again, and pull back the shell membrane and chorioallantoic membrane to expose the embryo.

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4-Resterilize the forceps, grasp the embryo loosely around the neck, and remove the entire embryo from the egg to a sterile Petri dish.

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5-Using forceps + scissors decapitate and eviscerate the embryo. Mince the embryo carcass into very small fragment with scissors .

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6- Add about 10ml sterile saline A to tissue fragment in the Petri dish,

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Swirl gently for 1 to 2 minutes, and carefully pour the entire content into a 125 ml Erlenmeyer flask ,tilt flask ,for decant saline A. Discard the saline A.

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7-Add 10 ml of sterile warm trypsin and pour into centrifuge filtrate at 1000 rpm for 10 minutes ,discard

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supernatant and resuspend the cells in 20 ml growth medium .

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8-Mix well for counting in a hemacytometer.

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9-Dilute 0.1ml of cell suspension with 0.9ml of MEM and count cell in a hemacytometer

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10-Adjust concentration to 5×10^5 cells per ml growth medium and dispense into cultures tubes.

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11-Incubate at 37°C until monolayered (6-7 days).

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Note: be sure examine cell culture each day ,if the medium is yellow may need to be adjust by $7.5\% \text{NaHCO}_3$.if floating cell are present or the color of the medium indicate a basic PH the medium should be changed.

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