# MEGBI Training Courses

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**Institute for Genetic Engineering, Ecology and Health (IGEEH)**

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**Middle East Genetics and Biotechnology Institute (MEGBI)**

Main Road, Ras Nhache, Batroun, Lebanon,

[www.aecenar.com/institutes/megbi](http://www.aecenar.com/institutes/megbi)

Email: [info@aecenar.com](mailto:info@aecenar.com)

Member Institutes of

**AECENAR**

Association for Economical and Technological Cooperation  
in the Euro-Asian and North-African Region

[www.aecenar.com](http://www.aecenar.com)

**Remark: All photos in the 4 parts are taken from laboratory work in MEGBI, Ras Nhache, Lebanon.**

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5 days training course for researchers and technical staff

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In this course some fundamental techniques for working with DNA are trained. The duration of this course is 5 days.

Molecular biological and genetical experimental strategies are used today for many tasks in agriculture, pharmaceutical and medical sciences. For this the goal of this course block is to teach basic techniques in this field.

The techniques learned in this course block will be used, insha Allah, to clone a human SRY gene, starting with the isolation of the DNA, with PCR and then compiling a plasmid vector. In one of the next courses, this vector shall be used to transform a plant.

Learning goals:

- Gene isolation
- PCR
- Methods to produce and analyse recombinant DNA (Electrophoresis, Plasmid technology, Restriction enzymes)

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##### List of kits/materials from Qiagen

Name	Details	Cat. No.	List Price (for Lebanon)
<b>QIAquick Gel Extraction Kit (50)</b>	50 QIAquick Spin Columns, Buffers, Collection Tubes (2 ml)	28704	More than 100 EUR
<b>QIAprep Spin Miniprep Kit (50)</b>	For 50 high-purity plasmid minipreps: 50 QIAprep Spin Columns, Reagents, Buffers, Collection Tubes (2 ml)	27104	More than 100 EUR

Alternatively there can be used similar kits for example from Peqlab.

#### **Promega**

##### Contact in Lebanon

Promega distributor in Lebanon:

**Atom Medical Company, s.a.r.l.** , Mr. Tonus Ajjub, 03 628732

P.O. Box 90-1203, Dora, St. Joseph Street, Zohour Building, Beirut,

Tel: (961) 1 249836, Fax: (961) 1 249838, E-mail Address: [ATOMMED@cyberia.net.lb](mailto:ATOMMED@cyberia.net.lb)

##### List of kits/materials from Promega

Ligation kit:

pGEM®-T Easy Vector System II  A1380 20 reactions

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**Molecular cloning** refers to the procedure of isolating a defined DNA sequence and obtaining multiple copies of it *in vivo*. Cloning is frequently employed to amplify DNA fragments containing genes, but it can be used to amplify any DNA sequence such as promoters, non-coding sequences, chemically synthesised oligonucleotides and randomly fragmented DNA. Cloning is utilized in a wide array of biological experiments and technological applications such as large scale protein production.

In the classical restriction and ligation cloning protocols, cloning of any DNA fragment essentially involves the following steps:

1. Isolation of the genetic material (DNA or cDNA) from the donor organism
2. DNA fragmentation with restriction endonucleases or specific amplification of DNA parts with PCR
3. ligation of DNA fragments to a plasmid (vector)
4. transformation

(transfection (Putting vector into a cell, which is therefor transformed), ...)  
screening/selection

Although these steps are invariable among cloning procedures a number of alternative routes can be selected at various points depending on the particular application; these are summarized as a 'cloning strategy'.

#### Isolation of insert (1. and 2.)

Initially, the DNA fragment to be cloned needs to be isolated. Preparation of DNA fragments for

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
cloning can be accomplished in a number of alternative ways. Insert preparation is frequently achieved by means of polymerase chain reaction, but it may also be accomplished by restriction enzyme digestion, DNA sonication and fractionation by agarose gel electrophoresis. Chemically synthesized oligonucleotides can also be used if the target sequence size does not exceed the limit of chemical synthesis. Isolation of insert can be done by using shotgun cloning, c-DNA clones, gene machines (artificial chemical synthesis).

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#### Ligation (3.)

#### Transformation (4.)

Following ligation, the ligation product (plasmid) is transformed into bacteria for propagation. The bacteria is then plated on selective agar to select for bacteria that has your plasmid of interest. Individual colonies are picked and tested for the wanted insert. Maxiprep can be done to obtain large quantity of the plasmid containing your inserted gene.

Transfection (Putting vector into a cell, which is therefor transformed)

Following ligation, a portion of the ligation reaction, including vector with insert in the desired orientation is transfected into cells. A number of alternative techniques are available, such as chemical sensitization of cells, electroporation and biolistics. Chemical sensitization of cells is frequently employed since this does not require specialized equipment and provides relatively high transformation efficiencies. Electroporation is used

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when extremely high transformation efficiencies are required, as in very inefficient cloning strategies. Biolistics are mainly utilized in plant cell transformations, where the cell wall is a major obstacle in DNA uptake by cells. The bacterial transformation is generally observed by blue white screening.

#### Selection (5.)

Finally, the transfected cells are cultured. As the aforementioned procedures are of particularly low efficiency, there is a need to identify the cells that contain the desired insert at the appropriate orientation and isolate these from those not successfully transformed. Modern cloning vectors include selectable markers (most frequently antibiotic resistance markers) that allow only cells in which the vector, but not necessarily the insert, has been transfected to grow. Additionally, the cloning vectors may contain colour selection markers which provide blue/white screening (via  $\alpha$ -factor complementation) on X-gal medium. Nevertheless, these selection steps do not absolutely guarantee that the DNA insert is present in the cells. Further investigation of the resulting colonies is required to confirm that cloning was successful. This may be accomplished by means of PCR, restriction fragment analysis and/or DNA sequencing.


Figure: Transformation and selection

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
**Deoxyribonucleic acid (DNA)** is a nucleic acid that contains the genetic instructions used in the development and functioning of all known living organisms and some viruses. The main role of DNA molecules is the long-term storage of information. DNA is often compared to a set of blueprints or a recipe, since it contains the instructions needed to construct other components of cells, such as proteins and RNA molecules. The DNA segments that carry this genetic information are called genes, but other DNA sequences have structural purposes, or are involved in regulating the use of this genetic information.

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Chemically, DNA consists of two long polymers of simple units called nucleotides, with backbones made of sugars and phosphate groups joined by ester bonds. These two strands run in opposite directions to each other and are therefore anti-parallel. Attached to each sugar is one of four types of molecules called bases. It is the sequence of these four bases along the backbone that encodes

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information. This information is read using the genetic code, which specifies the sequence of the amino acids within proteins. The code is read by copying stretches of DNA into the related nucleic acid RNA, in a process called transcription.

Within cells, DNA is organized into structures called chromosomes. These chromosomes are duplicated before cells divide, in a process called DNA replication. Eukaryotic organisms (animals, plants, fungi, and protists) store their DNA inside the cell nucleus, while in prokaryotes (bacteria and archae) it is found in the cell's cytoplasm. Within the chromosomes, chromatin proteins such as histones compact and organize DNA. These compact structures guide the interactions between DNA and other proteins, helping control which parts of the DNA are transcribed

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Organism	Genome size (base pairs)	Note
<u>Virus, Phage λ</u>	50,000	
<u>Bacterium, Escherichia coli</u>	4,000,000	
<u>Plant, Arabidopsis thaliana</u>	157,000,000	First plant genome sequenced, Dec 2000.
<u>Yeast, Saccharomyces cerevisiae</u>	20,000,000	

<u>Nematode, <i>Caenorhabditis elegans</i></u>	98,000,000	First multicellular animal genome, December 1998 <sup>[11]</sup>
<u>Insect, <i>Drosophila melanogaster</i></u> aka Fruit Fly	130,000,000	
<u>Fish, <i>Tetraodon nigroviridis</i></u> , type of Puffer fish	385,000,000	Smallest vertebrate genome known
Human	3,200,000,000	

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The polymerase chain reaction (PCR) is a technique widely used in molecular biology. It derives its name from one of its key components, a DNA polymerase used to amplify a piece of DNA by in vitro enzymatic replication. As PCR progresses, the DNA thus generated is itself used as a template for replication. This sets in motion a chain reaction in which the DNA template is exponentially amplified. With PCR it is possible to amplify a single or few copies of a piece of DNA across several orders of magnitude, generating millions or more copies of the DNA piece. PCR can be extensively modified to perform a wide array of genetic manipulations.

PCR is very versatile. Many types of samples can be analyzed for nucleic acids. Most PCR uses DNA as a target, rather than RNA, because of the stability of the DNA molecule and the ease with which DNA can be isolated. Almost all PCR applications employ a heat-stable DNA polymerase, such as Taq polymerase,

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an enzyme originally isolated from the bacterium *Thermus aquaticus*. This DNA polymerase enzymatically assembles a new DNA strand from DNA building blocks, the nucleotides, using single-stranded DNA as template and DNA oligonucleotides (also called DNA primers) required for initiation of DNA synthesis. The vast majority of PCR methods use thermal cycling, i.e., alternately heating and cooling the PCR sample to a defined series of temperature steps. These thermal cycling steps are necessary to physically separate the strands (at high temperatures) in a DNA double helix (DNA melting) used as template during DNA synthesis (at lower temperatures) by the DNA polymerase to selectively amplify the target DNA. The selectivity of PCR results from the use of primers that are complementary to the DNA region targeted for amplification under specific thermal cycling conditions.


Developed in 1983 by Kary Mullis, PCR is now a common and often indispensable technique used in medical and biological research labs for a variety of applications. These include DNA cloning for sequencing, DNA-based phylogeny, or functional analysis of genes; the diagnosis of hereditary diseases; the identification of genetic fingerprints (used in forensic sciences and paternity testing); and the detection and diagnosis of infectious diseases. In 1993 Mullis won the Nobel Prize in Chemistry for his work on PCR.

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PCR is used to amplify specific regions of a DNA strand (the DNA target). This can be a single gene, a part of a gene, or a non-coding sequence. Most PCR methods typically amplify DNA fragments of up to 10 kilo base pairs (kb), although some techniques allow for amplification of fragments up to 40 kb in size

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A basic PCR set up requires several components and reagents. These components include:

*DNA template* that contains the DNA region (target) to be amplified.

Two *primers*, which are complementary to the DNA regions at the 5' (five prime) or 3' (three prime) ends of the DNA region.

A DNA polymerase such as *Taq polymerase* or another DNA polymerase with a temperature optimum at around 70°C.

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*Deoxynucleoside triphosphates* (dNTPs; also very commonly and erroneously called deoxynucleotide triphosphates), the building blocks from which the DNA polymerases synthesizes a new DNA strand.

*Buffer solution*, providing a suitable chemical environment for optimum activity and stability of the DNA polymerase.

*Divalent* cations, magnesium or manganese ions; generally  $Mg^{2+}$  is used, but  $Mn^{2+}$  can be utilized for PCR-mediated DNA mutagenesis, as higher  $Mn^{2+}$  concentration increases the error rate during DNA synthesis.

*Monovalent cation* potassium ions.

The PCR is commonly carried out in a reaction volume of 10-200  $\mu$ l in small reaction tubes (0.2-0.5 ml volumes) in a thermal cycler. The thermal cycler heats and cools the reaction tubes to achieve the temperatures required at each step of the reaction (see below). Many modern thermal cyclers make use of the Peltier effect which permits both heating and cooling of the block holding the PCR tubes simply by reversing the electric current. Thin-walled reaction tubes permit favorable thermal conductivity to allow for rapid thermal equilibration. Most thermal cyclers have heated lids to prevent condensation at the top of the reaction tube. Older thermocyclers lacking a heated lid require a layer of oil on top of the reaction mixture or a ball of wax inside the tube.

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
**Figure 2:** Schematic drawing of the PCR cycle. (1) **Denaturing at 94-96°C.** (2) **Annealing at ~65°C** (3) **Elongation at 72°C.** Four cycles are shown here. The blue lines represent the DNA template to which primers (red arrows) anneal that are extended by the DNA polymerase (light green circles), to give shorter DNA products (green lines), which themselves are used as templates as PCR progresses. The PCR usually consists of a series of 20 to 40 repeated temperature changes called cycles; each cycle typically consists of 2-3 discrete temperature steps. Most commonly PCR is carried out with cycles that have three temperature steps (Fig. 2).

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The cycling is often preceded by a single temperature step (called *hold*) at a high temperature (>90°C), and followed by one hold at the end for final product extension or brief storage. The temperatures used and the length of time they are applied in each cycle depend on a variety of parameters. These include the enzyme used for DNA synthesis, the concentration of divalent ions and dNTPs in the reaction, and the melting temperature ( $T_m$ ) of the primers.

*Initialization step:* This step consists of heating the reaction to a temperature of 94-96°C (or 98°C if extremely thermostable

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polymerases are used), which is held for 1-9 minutes. It is only required for DNA polymerases that require heat activation by hot-start PCR.

*Denaturation step:* This step is the first regular cycling event and consists of heating the reaction to 94-98°C for 20-30 seconds. It causes melting of DNA template and primers by disrupting the hydrogen bonds between complementary bases of the DNA strands, yielding single strands of DNA.

*Annealing step:* The reaction temperature is lowered to 50-65°C for 20-40 seconds allowing annealing of the primers to the single-stranded DNA template. Typically the annealing temperature is about 3-5 degrees Celsius below the T<sub>m</sub> of the primers used. Stable DNA-DNA hydrogen bonds are only formed when the primer sequence very closely matches the template sequence. The polymerase binds to the primer-template hybrid and begins DNA synthesis.

*Extension/elongation step:* The temperature at this step depends on the DNA polymerase used; Taq polymerase has its optimum activity temperature at 75-80°C,<sup>[10][11]</sup> and commonly a temperature of 72°C is used with this enzyme. At this step the DNA polymerase synthesizes a new DNA strand complementary to the DNA template strand by adding dNTPs that are complementary to the template in 5' to 3' direction, condensing the 5'-phosphate group of the dNTPs with the 3'-hydroxyl group at the end of the nascent (extending) DNA strand. The extension time depends both on the DNA polymerase used and on the length of the DNA fragment to be amplified. As a rule-of-thumb, at its optimum temperature, the DNA polymerase will polymerize a thousand bases per minute. Under optimum conditions, i.e., if there are no limitations due to limiting substrates or reagents, at each extension step, the amount of DNA target is doubled, leading to exponential (geometric) amplification of the specific

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DNA fragment.

*Final elongation:* This single step is occasionally performed at a temperature of 70-74°C for 5-15 minutes after the last PCR cycle to ensure that any remaining single-stranded DNA is fully extended.


*Final hold:* This step at 4-15°C for an indefinite time may be employed for short-term storage of the reaction.

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**Figure 3:** Ethidium bromide-stained PCR products after gel electrophoresis. Two sets of primers were used to amplify a target sequence from three different tissue samples. No amplification is present in sample #1; DNA bands in sample #2 and #3 indicate successful amplification of the target sequence. The gel also shows a positive control, and a DNA ladder containing DNA fragments of defined length for sizing the bands in the experimental PCRs.

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To check whether the PCR generated the anticipated DNA fragment (also sometimes referred to as the amplicon or amplicon), agarose gel

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electrophoresis is employed for size separation of the PCR products. The size(s) of PCR products is determined by comparison with a DNA ladder (a molecular weight marker), which contains DNA fragments of known size, run on the gel alongside the PCR products (see Fig. 3).

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The PCR process can be divided into three stages:

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*Exponential amplification:* At every cycle, the amount of product is doubled (assuming 100% reaction efficiency). The reaction is very specific and precise.

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*Levelling off stage:* The reaction slows as the DNA polymerase loses activity and as consumption of reagents such as dNTPs and primers causes them to become limiting.

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*Plateau:* No more product accumulates due to exhaustion of reagents and enzyme.

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In practice, PCR can fail for various reasons, in part due to its sensitivity to contamination causing amplification of spurious DNA products. Because of this, a number of techniques and procedures have been developed for optimizing PCR conditions. Contamination with extraneous DNA is addressed with lab protocols and procedures that separate pre-PCR mixtures from potential DNA contaminants. This usually involves spatial separation of PCR-setup areas from areas for analysis or purification of PCR products, and thoroughly cleaning the work surface between reaction setups. Primer-design techniques are important in improving PCR product yield and in avoiding the formation of spurious products, and the usage of alternate buffer components or polymerase enzymes can help with amplification of

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PCR allows isolation of DNA fragments from genomic DNA by selective amplification of a specific region of DNA. This use of PCR augments many methods, such as generating hybridization probes for Southern or northern hybridization and DNA cloning, which require larger amounts of DNA, representing a specific DNA region. PCR supplies these techniques with high amounts of pure DNA, enabling analysis of DNA samples even from very small amounts of starting material.


Other applications of PCR include DNA sequencing to determine unknown PCR-amplified sequences in which one of the amplification primers may be used in Sanger sequencing, isolation of a DNA sequence to expedite recombinant DNA technologies involving the insertion of a DNA sequence into a plasmid or the genetic material of another organism. Bacterial colonies (E.coli) can be rapidly screened by PCR for correct DNA vector constructs.

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**Figure 4:** Electrophoresis of PCR-amplified DNA fragments. (1) Father. (2) Child. (3) Mother. The child has inherited some, but not all of the fingerprint of each of its parents, giving it a new, unique fingerprint.

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PCR may also be used for *genetic fingerprinting*; a forensic technique used to identify a person or organism by comparing experimental DNAs through different PCR-based methods. Some PCR 'fingerprints' methods have high discriminative power and can be used to identify genetic relationships between individuals, such as parent-child or between siblings, and are used in paternity testing (Fig. 4). This technique may also be used to determine evolutionary relationships among organisms.

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Because PCR amplifies the regions of DNA that it targets, PCR can be used to analyze extremely small amounts of sample. This is often critical for forensic analysis, when only a trace amount of DNA is available as evidence. PCR may also be used in the analysis of ancient DNA that is tens of thousands of years old. These PCR-based techniques have been successfully used on animals, such as a forty-thousand-year-old mammoth, and also on human DNA, in applications ranging from the analysis of Egyptian mummies to the identification of a Russian Tsar.

Quantitative PCR methods allow the estimation of the amount of a given sequence present in a sample – a technique often applied to quantitatively determine levels of gene expression. Real-time PCR is an established tool for DNA quantification that measures the accumulation of DNA product after each round of PCR amplification.

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PCR allows early diagnosis of malignant diseases such as leukemia and lymphomas, which is currently the highest developed in cancer research and is already being used routinely. PCR assays can be performed directly on genomic DNA samples to detect translocation-specific malignant cells at a sensitivity which is at least 10,000 fold higher than other methods.

PCR also permits identification of non-cultivable or slow-growing microorganisms such as mycobacteria, anaerobic bacteria, or viruses from tissue culture assays and animal models. The basis for PCR diagnostic applications in microbiology is the detection of infectious agents and the discrimination of non-pathogenic from pathogenic strains by virtue of specific genes.

Viral DNA can likewise be detected by PCR. The primers used need to be specific to the targeted sequences in the DNA of a virus, and the PCR can be used for diagnostic analyses or DNA sequencing of the viral genome. The high sensitivity of PCR permits virus detection soon after infection and even before the onset of disease. Such early detection may give physicians a significant lead in treatment. The amount of virus ("viral load") in a patient can also be quantified by PCR-based DNA quantitation techniques (see below).

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A 1971 paper in the Journal of Molecular Biology by Kleppe and co-workers first described a method using an enzymatic assay to replicate a short DNA template with primers *in vitro*. However, this early manifestation of the basic PCR principle did not receive much attention, and the invention of the polymerase chain reaction in 1983 is generally credited to Kary Mullis. At the core of the PCR method is the use of

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a suitable DNA polymerase able to withstand the high temperatures of  $>90^{\circ}\text{C}$  ( $>195^{\circ}\text{F}$ ) required for separation of the two DNA strands in the DNA double helix after each replication cycle. The DNA polymerases initially employed for in vitro experiments presaging PCR were unable to withstand these high temperatures. So the early procedures for DNA replication were very inefficient, time consuming, and required large amounts of DNA polymerase and continual handling throughout the process.

A 1976 discovery of Taq polymerase a DNA polymerase purified from the thermophilic bacterium, *Thermus aquaticus*, which naturally occurs in hot ( $50$  to  $80^{\circ}\text{C}$  ( $120$  to  $175^{\circ}\text{F}$ )) environments paved the way for dramatic improvements of the PCR method. The DNA polymerase isolated from *T. aquaticus* is stable at high temperatures remaining active even after DNA denaturation, thus obviating the need to add new DNA polymerase after each cycle. This allowed an automated thermocycler-based process for DNA amplification.

At the time he developed PCR in 1983, Mullis was working in Emeryville, California for Cetus Corporation, one of the first biotechnology companies. There, he was responsible for synthesizing short chains of DNA. Mullis has written that he conceived of PCR while cruising along the Pacific Coast Highway one night in his car. He was playing in his mind with a new way of analyzing changes (mutations) in DNA when he realized that he had instead invented a method of amplifying any DNA region through repeated cycles of duplication driven by DNA polymerase.

In *Scientific American*, Mullis summarized the procedure: "Beginning with a single molecule of the genetic material DNA, the PCR can generate 100 billion similar molecules in an afternoon. The reaction is easy to execute. It requires no more than a test tube, a few simple reagents, and a source of heat." He was awarded the Nobel

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Prize in Chemistry in 1993 for his invention, seven years after he and his colleagues at Cetus first put his proposal to practice. However, some controversies have remained about the intellectual and practical contributions of other scientists to Mullis' work, and whether he had been the sole inventor of the PCR principle. (see main article: Kary Mullis)

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The PCR technique was patented by Cetus Corporation, where Mullis worked when he invented the technique in 1983. The *Taq* polymerase enzyme was also covered by patents. There have been several high-profile lawsuits related to the technique, including an unsuccessful lawsuit brought by DuPont. The pharmaceutical company Hoffmann-La Roche purchased the rights to the patents in 1992 and currently holds those that are still protected.

A related patent battle over the *Taq* polymerase enzyme is still ongoing in several jurisdictions around the world between Roche and Promega. The legal arguments have extended beyond the life of the original PCR and *Taq* polymerase patents, which expired on March 28, 2005

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A **restriction enzyme** (or **restriction endonuclease**) is an enzyme that cuts double-stranded DNA following its specific recognition of short nucleotide sequences, known as restriction sites, in the DNA. Such enzymes, found in bacteria and archaea, are thought to have evolved to provide a defense mechanism against invading viruses. Inside a bacterial host, the restriction enzymes selectively cut up *foreign* DNA in a process called *restriction*; host DNA is methylated by a modification enzyme (a methylase) to protect it from the restriction enzyme's activity. Collectively, these two processes form the restriction modification system. To cut the DNA, a restriction enzyme makes two incisions, once through each sugar-phosphate backbone (i.e. each strand) of the DNA double helix.

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### Restriction enzyme glossary

- **restriction** – is a specific type of rule that defines a finite (and generally absolute) boundary defined for a type of process or function
- **enzyme** – a protein that catalyzes a chemical reaction
- **endonuclease** – an enzyme that cleaves the phosphodiester bond within a polynucleotide chain such as DNA. A restriction enzyme is a type of endonuclease.
- **molecular recognition** – used by restriction enzymes to locate specific sequences of DNA on which to bind and subsequently cleave
- **recognition site** or **recognition sequence** – the DNA location/sequence to which restriction enzymes bind
- **restriction site** – the DNA sequence that is cleaved by the restriction enzyme

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5'-GTATAC-3'           3'-CATATG-5'	A palindromic recognition site reads the same on the reverse strand as it does on the forward strand
---	--

Restriction enzymes recognize a specific sequence of nucleotides<sup>[2]</sup> and produce a double-stranded cut in the DNA that prevents the phage from replicating. While recognition sequences vary widely, with lengths between 4 and 8 nucleotides, many of them are palindromic; that is, the sequence on one strand reads the same in the reverse direction on the complementary strand.<sup>[13]</sup> The meaning of "palindromic" in this context is different from what one might expect from its linguistic usage: GTAATG is not a palindromic DNA sequence, but GTATAC is (GTATAC is complementary to CATATG).

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Recognition sequences in DNA differ for each restriction enzyme, producing differences in the length, sequence and strand orientation (5' end or the 3' end) of a sticky-end "overhang" of an enzyme restriction. Different restriction enzymes that

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recognize the same sequence are known as neoschizomers. These often cleave in a different location of the sequence; however, the specific type that cleaves in the same location as the prototype is known as an isoschizomer.

Bacteria prevent their own DNA from being cut by modifying their nucleotides via DNA methylation.

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Restriction endonucleases are categorized into three general groups (Types I, II and III) based on their composition and enzyme cofactor requirements, the nature of their target sequence, and the position of their DNA cleavage site relative to the target sequence.

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Type I restriction enzymes were the first to be identified and are characteristic of two different strains (K-12 and B) of *E. coli*.<sup>[18]</sup> These enzymes cut at a site that differs, and is some distance (at least 1000 bp) away, from their recognition site. The recognition site is asymmetrical and is composed of two portions – one containing 3-4 nucleotides, and another containing 4-5 nucleotides – separated by a spacer of about 6-8 nucleotides. Several enzyme cofactors, including S-Adenosyl methionine (AdoMet), hydrolyzed adenosine triphosphate (ATP) and magnesium ( $Mg^{2+}$ ) ions, are required for their activity. Type I restriction enzymes possess three subunits called HsdR, HsdM, and HsdS; HsdR is required for restriction, HsdM is necessary for adding methyl groups to host DNA (methyltransferase activity) and HsdS is important for specificity of cut site recognition in addition to its methyltransferase activity

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Typical type II restriction enzymes differ from type I restriction enzymes in several ways. They are composed of only one subunit, their recognition sites are usually undivided and palindromic and 4-8 nucleotides in length, they recognize and cleave DNA at the same site, and they do not use ATP or AdoMet for their activity – they usually require only  $Mg^{2+}$  as a cofactor. These are the most commonly available and used restriction enzymes. In the 1990s and early 2000s, new enzymes from this family were discovered that did not follow all the classical criteria of this enzyme class, and new subfamily nomenclature was developed to divide this large family into subcategories based on deviations from typical characteristics of type II enzymes. These subgroups are defined using a letter suffix.

Type IIB restriction enzymes (e.g. Bcgl and BpII) are multimers, containing more than one subunit. They cleave DNA on both sides of their recognition to cut out the recognition site. They require both AdoMet and  $Mg^{2+}$  cofactors. Type IIE restriction endonucleases (e.g. NaeI) cleave DNA following interaction with two copies of their recognition sequence. One recognition site acts as the target for cleavage, while the other acts as an allosteric effector that speeds up or improves the efficiency of enzyme cleavage. Similar to type IIE enzymes, type IIF restriction endonucleases (e.g. NgoMIV) interact with two copies of their recognition sequence but cleave both sequences at the same time. Type IIG restriction endonucleases (Eco57I) do have a single subunit, like classical Type II restriction enzymes, but require the cofactor AdoMet to be active. Type IIM restriction endonucleases, such as DpnI, are able to recognize and cut methylated DNA. Type IIS restriction endonucleases (e.g. FokI) cleave DNA at a defined distance from their non-palindromic asymmetric recognition sites. These enzymes may function as dimers.

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Similarly, Type III restriction enzymes (e.g., Bpu10I and BslI) are composed of two different subunits. Some recognize palindromic sequences while others have asymmetric recognition sites.

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Type III restriction enzymes (e.g. *EcoP15*) recognize two separate non-palindromic sequences that are inversely oriented. They cut DNA about 20-30 base pairs after the recognition site.<sup>[19]</sup> These enzymes contain more than one subunit and require AdoMet and ATP cofactors for their roles in DNA methylation and restriction, respectively.

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Since their discovery in the 1970s, more than 100 different restriction enzymes have been identified in different bacteria. Each enzyme is named after the bacterium from which it was isolated using a naming system based on bacterial genus, species and strain. For example, the name of the *EcoRI* restriction enzyme was derived as shown in the box.

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I First (order of identification in identified the bacterium)

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Isolated restriction enzymes are used to manipulate DNA for different scientific applications.

They are used to assist insertion of genes into plasmid vectors during gene cloning and protein expression experiments. For optimal use, plasmids that are commonly used for gene cloning are modified to include a short *linker* sequence (called the multiple cloning site,

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or MCS) rich in restriction enzyme recognition sequences. This allows flexibility when inserting gene fragments into the plasmid vector; restriction sites contained naturally within genes influence the choice of endonuclease for digesting the DNA since it is necessary to avoid restriction of wanted DNA while intentionally cutting the ends of the DNA. To clone a gene fragment into a vector, both plasmid DNA and gene insert are typically cut with the same restriction enzymes, and then glued together with the assistance of an enzyme known as a DNA ligase.

They can be used to distinguish gene alleles by specifically recognizing single base changes in DNA known as single nucleotide polymorphisms (SNPs). This is only possible if a SNP alters the restriction site present in the allele. In this method, the restriction enzyme can be used to genotype a DNA sample without the need for expensive gene sequencing. The sample is first digested with the restriction enzyme to generate DNA fragments, and then the different sized fragments separated by gel electrophoresis. In general, alleles with correct restriction sites will generate two visible bands of DNA on the gel, and those with altered restriction sites will not be cut and will generate only a single band. The number of bands reveals the sample subject's genotype, an example of restriction mapping.

In a similar manner, restriction enzymes are used to digest genomic DNA for gene analysis by Southern Blot. This technique allows researchers to identify how many copies (or paralogues) of a gene are present in the genome of one individual, or how many gene mutations (polymorphisms) have occurred within a population. The latter example is called Restriction Fragment Length Polymorphism (RFLP).

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Examples of restriction enzymes include:

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
Figure: The pGEX-3x plasmid is a popular cloning vector.

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A **cloning vector** is a small DNA vehicle into which a foreign DNA fragment can be inserted. The insertion of the fragment into the cloning vector is carried out by treating the vehicle and the foreign DNA with the same restriction enzyme, then ligating the fragments together. There are many types of cloning vectors. Genetically engineered plasmids and bacteriophages (such as phage  $\lambda$ ) are perhaps most commonly used for this purpose. Other types of cloning vectors include bacterial artificial chromosomes (BACs) and yeast artificial chromosomes (YACs).

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Most commercial cloning vectors have key features that have made their use in molecular biology so widespread.

In the case of expression vectors, the main purpose of these vehicles is the controlled expression of a particular gene inside a convenient host organism (eg. *E. coli*). Control of expression can be very important; it is usually desirable to insert the target DNA into a site that is under the control of a particular promoter. Some commonly used promoters are T7 promoters, *lac* promoters (*bla* promoter) and cauliflower mosaic virus's 35s promoter (for plant vectors).

To allow for convenient and favorable insertions, most cloning vectors have had nearly all their restriction sites engineered

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out of them and a synthetic multiple cloning site (MCS) inserted that contains many restriction sites. MCSs allow for insertions of DNA into the vector to be targeted and possibly directed in a chosen orientation. A selectable marker, such as an antibiotic resistance [eg. beta-lactamase (see figure)] is often carried by the vector to allow the selection of positively transformed cells (see Screening below). All plasmids must carry a functional origin of replication (ORI; not shown in figure). Some other possible features present in cloning vectors are: *vir* genes for plant transformation, intergrase sites for chromosomal insertion, *lacZα* fragment for  $\alpha$  complementation and blue-white selection, and/or reporter genes in frame with and flanking the MCS to facilitate the production of recombinant proteins [eg. fused to the Green fluorescent protein (GFP) or to the glutathione S-transferase (see figure)].

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Many general purpose vectors such as pUC19 usually include a system for detecting the presence of a cloned DNA fragment, based on the loss of an easily scored phenotype. The most widely used is the gene coding for *E. coli*  $\beta$ -galactosidase, whose integrity can easily be detected by the ability of the enzyme it encodes to hydrolyze the soluble, colourless substrate X-gal (5-bromo-4-chloro-3-indolyl-beta-d-galactoside) into an insoluble, blue product (5,5'-dibromo-4,4'-dichloro indigo). Cloning a fragment of DNA within the vector-based gene encoding the  $\beta$ -galactosidase prevents the production of an active enzyme. If X-gal is included in the selective agar plates, transformant colonies are generally blue in the case of a vector with no inserted DNA and white in the case of a vector containing a fragment of cloned DNA.

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In molecular biology, **DNA ligase** is a particular type of ligase (EC 6.5.1.1) that can link together two DNA strands that have single-strand breaks (a break in both complementary strands of DNA). The alternative, a double-strand break, is fixed by a different type of DNA ligase using the complementary strand as a template but still requires DNA ligase to create the final phosphodiester bond to fully repair the DNA.

DNA ligase has applications in both DNA repair and DNA replication. In addition, DNA ligase has extensive use in molecular biology laboratories for Genetic recombination experiments.

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
The mechanism of DNA ligase is to form two covalent phosphodiester bonds between 3' hydroxyl ends of one nucleotide with the 5' phosphate end of another. ATP is not required for the ligase reaction. A pictorial example of how a ligase works (with sticky ends):

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Ligase will also work with blunt ends, although higher enzyme concentrations and different reaction conditions are required.

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DNA ligases have become an indispensable tool in modern molecular biology research for generating recombinant DNA sequences. For example, DNA ligases are used with restriction enzymes to insert DNA fragments, often genes, into plasmids. One vital, and often tricky, aspect to performing successful recombination experiments involving ligase is controlling the optimal temperature. Most experiments use T4 DNA Ligase (isolated from bacteriophage T4) which is most active at 25°C. However in order to perform successful ligations, the optimal enzyme temperature needs to be balanced with the melting temperature  $T_m$  (also the annealing temperature) of the DNA fragments being ligated.

If the ambient temperature exceeds  $T_m$ , homologous pairing of the sticky ends will not occur because the high temperature disrupts hydrogen bonding. The shorter the DNA fragments, the lower the  $T_m$ . Thus for sticky ends (overlaps) less than ten base pairs long, ligation experiments are performed at very low temperatures (~4-8°C) for a long period of time (often overnight).

The common commercially available DNA ligases were originally discovered in bacteriophage T4, *E. coli* and other bacteria.

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DNA ligation is the process of joining together two DNA molecule ends (either from the same or different molecules). Specifically, it involves creating a phosphodiester bond between the 3' hydroxyl of one nucleotide and the 5' phosphate of another. This reaction is usually catalyzed by a DNA ligase enzyme. This enzyme will ligate DNA fragments having blunt or overhanging, complementary, 'sticky' ends. Typically, it is easier to ligate molecules with

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complementary sticky ends than blunt ends. T4 DNA ligase is the most commonly used DNA ligase for molecular biology techniques and can ligate 'sticky' or blunt ends.

The two components of the DNA in the ligation reaction should be equimolar and around 100µg/ml. Most commonly, one wants to ligate an insert DNA molecule into a plasmid, ready for bacterial transformation. Typically, DNA and plasmid vector are individually cut to yield complementary ends, then both are added to a ligation reaction to be circularised by DNA ligase. If the plasmid backbone to insert DNA ratio is too high then excess 'empty' mono and polymeric plasmids will be generated. If the ratio is too low then the result may be an excess of linear and circular homo- and heteropolymers.

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$$\text{Insert Mass in ng} = 6 \times \left[ \frac{\text{Insert Length in bp}}{\text{Vector Length in bp}} \right] \times \text{Vector Mass in ng}$$

The insert to vector molar ratio can have a significant effect on the outcome of a ligation and subsequent transformation step. Molar ratios can vary from a 1:1 insert to vector molar ratio to 10:1. It may be necessary to try several ratios in parallel for best results.

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1. Add appropriate amount of deionized H<sub>2</sub>O to sterile 0.6 mL tube
2. Add 1 µL ligation buffer to the tube.  
Vortex buffer before pipetting to ensure that it is well-mixed. Remember that the buffer contains ATP so repeated freeze, thaw cycles can degrade the ATP thereby decreasing the efficiency of ligation.
3. Add appropriate amount of insert to the tube.
4. Add appropriate amount of vector to the tube.
5. Add 0.5 µL ligase. Vortex ligase before pipetting to ensure that it is well-mixed. Also, the ligase, like most enzymes, is in some percentage of glycerol which tends to stick to the sides of your tip. To ensure you add only 0.5 µL, just touch your tip to the surface of the liquid when pipetting.
6. Let the 10 µL solution sit at 22.5°C for 30 mins
7. Denature the ligase at 65°C for 10min
8. Dialyze for 20 minutes if electroporating
9. Use disks shiny side up

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A protocol analysis experiment for a typical DNA ligation (7.2 kb vector + 0.6 kb insert, sticky ends) gave optimal ligation efficiency when 50 ng of vector was ligated overnight at 16°C with a 2:1 insert:vector molar ratio and standard T4 ligase. Ligase was heat inactivated at 65°C for 20 mins before 2 µL (of 20 µL) was used to transform commercial heat-shock competent cells. Ligation efficiency was **marginally decreased** by doing a 1 hr ligation at room temperature using 100 ng vector

1. Using insert:vector molar ratios of 5:1 and 1:1

Ligation efficiency was **noticeably decreased** (x100) by

1. Sticky end ligation with a larger insert (5.2 kb vector + 2.6 kb insert)

2. Blunt end ligation

Ligation efficiency was **severely decreased** (x10000) by

1. Using DNA fragments that have been exposed to ethidium bromide and UV during the gel extraction procedure (*difficult to avoid but heartily recommended*)

2. Using the NEB Quick Ligation Kit (*possibly a bad batch*)

For additional trouble shooting, check out the NEB FAQ page for T4 ligation:

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1. Make sure the buffer is completely melted and dissolved. The white precipitate is BSA according to NEB. Make sure the buffer still smells strongly

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like "wet dog" (to check if the DTT is still good).

2. Because ligase buffer contains ATP, which is unstable and degraded by multiple freeze/thaw cycles, you may want to make 10-20ul aliquots from the original tube. Ligase buffer may be spiked with additional ATP.

3. If you are having trouble with your ligation, NEB offers FAQ's (Quick Ligation T4 DNA ligase) and tips to help.

4. Prior to the ligation, some heat their DNA slightly (maybe ~37°C) to melt any sticky ends which may have annealed improperly at low temperatures.

5. Tom Knight has read that ligase can inhibit transformation. By heat-inactivating the ligase, this inhibition can be avoided. However, according to the NEB FAQ, heat-inactivation of PEG (which is present in the ligation reaction) also inhibits transformation, therefore a spin-column purification is recommended prior to transformation if you are having problems.

6. Treating PCR products with proteinase K prior to restriction digest dramatically improves the efficiency of subsequent ligation reactions.

7. Using SYBR Safe DNA Gel Stain is a safer, non-carcinogenic alternative to ethidium bromide.

8. T4 DNA Ligase is very sensitive to shear, so spinning your ligation mix or vortexing it to mix it can affect your yields. Instead try mixing with the pipette tip or slowly resuspending the solution.

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This protocol is primarily based on Endy:DNA ligation using T4 DNA ligase.

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
Crowe JS, Cooper HJ, Smith MA, Sims MJ, Parker D, and Gewert D. *Improved cloning efficiency of polymerase chain reaction (PCR) products after proteinase K digestion.* Nucleic Acids Res 1991 Jan 11; 19(1) 184. PMID:2011503. [PubMed](#) [HubMed](#) [Crowe-NAR-1991]

Olivera BM and Lehman IR. *Linkage of polynucleotides through phosphodiester bonds by an enzyme from Escherichia coli*. Proc Natl Acad Sci U S A 1967 May; 57(5) 1426-33. pmid:5341238. PubMed HubMed [Olivera-PNAS-1967]

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Minipreparation of plasmid DNA and Ding dong DNA is a rapid, small-scale isolation of plasmid DNA from bacteria. It is based on the alkaline lysis method invented by the researchers Birnboim and Doly in 1979. The extracted plasmid DNA resulting from performing a miniprep is itself often called a "miniprep".

When bacteria are lysed under alkaline conditions both DNA and proteins are precipitated. Some scientists reduce the concentration of NaOH used to 0.1M in order to reduce the occurrence of ssDNA. After the addition of acetate-containing neutralization buffer the large and less supercoiled chromosomal DNA and proteins precipitate, but the small bacterial DNA plasmids can renature and stay in solution. Addition of phenol/chloroform can dissolve and denature proteins, like DNase. This is especially important if the plasmids are to be used for enzyme digestion. Otherwise, smearing may occur in enzyme restricted form of plasmid DNA.

Minipreps are used in the process of molecular cloning to analyze bacterial clones. A typical plasmid DNA yield of a miniprep is 20 to 30 µg depending on the cell strain.

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Typically 10-30 µl/sample of the DNA fragments to separate are obtained, as well as a mixture of DNA fragments (usually 10-20) of known size (after processing with

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DNA size markers either from a commercial source or prepared manually).

- Buffer solution, usually TBE buffer or TAE 1.0x, pH 8.0

- Agarose

- An ultraviolet-fluorescent dye, ethidium bromide, (5.25 mg/ml in H<sub>2</sub>O). The stock solution be careful handling this. Alternative dyes may be used, such as SYBR Green.

- zNitrile rubber gloves. Latex gloves do not protect well from ethidium bromide

- A color marker dye containing a low molecular weight dye such as "bromophenol blue" (to enable tracking the progress of the electrophoresis) and glycerol (to make the DNA solution more dense so it will sink into the wells of the gel).

- A gel rack

- A "comb"

- Power Supply

- UV lamp or UV lightbox or other method to visualize DNA in the gel

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There are several methods for preparing gels. A common example is shown here. Other methods might differ in the buffering system used, the sample size to be loaded, the total volume of the gel (typically thickness is kept to a constant amount while length and breadth are varied as needed). Most agarose gels used in modern biochemistry and molecular biology are prepared and run horizontally.

1 Make a 1% agarose solution in 100ml TAE, for typical DNA fragments (see figures). A solution of up to 2-4% can be used if you analyze small DNA molecules, and for large

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After the gel has been prepared, use a micropipette to inject about 2.5 µl of stained DNA (a DNA ladder is also highly recommended). Close the lid of the electrophoresis chamber and apply current (typically 100 V for 30 minutes with 15 ml of gel). The colored dye in the DNA ladder and DNA samples acts as a "front wave" that runs

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faster than the DNA itself. When the "front wave" approaches the end of the gel, the current is stopped. The DNA is stained with ethidium bromide, and is then visible under ultraviolet light.

1. The agarose gel with three slots/wells (S).
2. Injection of DNA ladder (**molecular weight** markers) into the first slot.
3. DNA ladder injected. Injection of samples into the second and third slot.
4. A current is applied. The DNA moves toward the positive anode due to the negative charges on its phosphate backbone.
5. Small DNA strands move fast, large DNA strands move slowly through the gel. The DNA is not normally visible during this process, so the marker dye is added to the DNA to avoid the DNA being run entirely off the gel. The marker dye has a low molecular weight, and migrates faster than the DNA, so as long as the marker has not run past the end of the gel, the DNA will still be in the gel.
6. Add the color marker dye to the DNA ladder.

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
Agarose gel with samples loaded in the slots, before the electrophoresis process

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A pattern of DNA-bands under UV light

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
Figure: Schematic drawing of the electrophoresis process, see text for description of steps

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#### 2.1.1 SRY gene / الجين SRY

SRY (Sex-determining Region Y) is a sex-determining gene on the Y chromosome in the therians (placental mammals and marsupials)

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This intronless gene encodes a transcription factor that is a member of the high mobility group (HMGtt)-box family of DNA-binding proteins. This protein is the testis determining factor (TDF), also referred to as the SRY protein, which initiates male sex determination.

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
The sequence is given in the direction 3'-----5'

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3D structure of TDF (there are several variants in the database, e.g. 1J46 and 1HRY):

Testis  
determining  
factor  
**TDF (HMG-  
Domain, green)  
with DNA (blue)  
according to  
PDB 1J46**

<http://www.rcsb.org/pdb/explore/images.do?structureId=1HRY>

FASTA format sequence of the protein (version 1hry):

```
>1HRY:A|PDBID|CHAIN|SEQUENCE
VQDRVKRP MNAFIVWSRDQRRKMALENPRMRNSEISKQLGYQWKMLTEAEKWPFFQEAQKLQAMHR
EKYPNYKYP
>1HRY:C|PDBID|CHAIN|SEQUENCE
GTTTGTGC
>1HRY:B|PDBID|CHAIN|SEQUENCE
GCACAAAC
```

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PBND buffer (for expose the membrane):

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Primers / hSRY5: 5'-gttgagggg tggtagggc ggaga-3'

المشروعات hSRY3: 5'-tgaaatg aataaggcct ttattagcca-3'

(hSRY3 is complementary sequence of 5'-tggttaataa aggcttatt catttca -3')

**Important remark:** The SRY gene has to be put later in pGEM-T vector and then has to be cut out from it afterwards with EcoRI. For this reason there has to be put at both primers a EcoRI cutting sequence. In this course primers without these EcoRI cutting sequences are used. So the last step of the course (restriction digest with EcoRI) is not possible.

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1. If you are male: Take a smear of your oral mucosa. Take it with a plastic spatula. Ask for a probe of a female oral mucosa smear.

If you are female: Take a smear of your oral mucosa. Take it with a plastic

Figure: Putting the smear probe in an Eppendorf tube

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spatula. Ask for a probe of a male oral mucosa smear

Each of the smear (marked as "male" and as "female") is flushed with 20µl PBNB buffer into an Eppendorf tube.

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For each Eppendorf tube:

2. Add 3µl Proteinase K (10mg/ml). Then incubate several (about 3-5) hours at 55 °C (for activation the enzyme).
3. Incubate at 95°C for 15 minutes (95°C: for denaturation the proteinase k).

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5. Performing the PCR reaction: Program the cyclor as follows:

40 cycles of:

40sec. at 95°C

40sec. at 53°C

1 min. at 72°C

(95°C for denaturation the double brind of DNA, 53°C for activation the primer, 72°C for activation the polymerase).

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1. TAE buffer 1.0x pH8.0 (prepared in part 2.1.5.4).

2. Agarose

3. Ethidium bromide (Stock solution: 1g Eth.bromide in 100 ml dest. water)

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- weight 1% agarose for 100ml TAE (see part 2.1.5.4).
- Dissolve the agarose on the agitator machine (agarose dissolve at 70 °C).
- Let the solution cool down to about 60°C at room temperature.
- Stir or swirl the solution while cooling.
- Put the gel on the gel rack and insert the comb.
- Let the gel to solidify.

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1. 20ml EDTA (0.5M).
2. 240ml acetic acid (5%).
3. 400 ml Tris base (0.01M).
4. 240ml distilled water.
5. Regulate with NaOH(2M) the pH on 8.2-8.4.

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
Figure: Gel rack with voltage supply

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I Take of each of your 3 Eppendorf tubes (male DNA, female DNA, DNA ladder) 4µl DNA and add as into a new Eppendorf tube as follows:

- For male smear probe: 4µl DNA (from smear) +3µl bromophenol blue + 3µl glycerol.
- For female smear probe: 4µl DNA (from smear) +3µl bromophenol blue + 3µl glycerol.
- For DNA ladder: 4µl ladder DNA +3µl bromophenol blue + 3µl glycerol.

II Inject 10µl from each eppendorf tube in the slots/wells of gel prepared on the gel rack

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III Close the lid of the electrophoresis chamber and apply the voltage of 80 V for 1h30min.

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IV After migration of the band; put the gel into a deep vessel (to put Ethidium bromide on it and to wash it afterwards).

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V Soaking of gel with Ethidium bromide: Put 250 ml dest. water and afterwards 30µl of ethidium bromide stock solution into the vessel with the gel. Let it soaking for 30 min.

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Warning: Be careful with ethidium bromide!!! Use specific gloves, put the water with ethidium bromide in a specific place, do not throw it in the nature. Ethidium bromide could cause cancer (propably onkogen) and it is mutagen.

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VI Wash the gel from ethidium bromide 3 times with dist. water (250ml).

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VII Put the gel on UV- machine. Insha Allah you will see bands for the DNA ladder and the male DNA:

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QIAquick Gel Extraction Kit Protocol using a microcentrifuge

This protocol is designed to extract and purify DNA of 70 bp to 10 kb from standard or low-melt agarose

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gels in TAE or TBE buffer. Up to 400 mg agarose can be processed per spin column. This kit can also be used for DNA cleanup from enzymatic reactions.

For DNA cleanup from enzymatic reactions using this protocol, add 3 volumes of Buffer QG and 1 volume of isopropanol to the reaction, mix, and proceed with step 6 of the protocol. Alternatively, use the MinElute Reaction Cleanup Kit.

#### Important points before starting

- The yellow color of Buffer QG indicates a pH  $\leq 7.5$
- Add ethanol (96–100%) to Buffer PE before use (see bottle label for volume).
- All centrifugation steps are carried out at  $17,900 \times g$  (13,000 rpm) in a conventional table-top microcentrifuge at room temperature.

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1. Cut the DNA fragment from the agarose gel with a clean sharp scalpel.

Minimize the size of the gel slice by removing extra agarose.

2. Weigh the gel slice in a colorless tube. Add 3 volumes of Buffer QG to 1 volume of gel (100 mg ~ 100  $\mu$ l).

For example, add 300  $\mu$ l of Buffer QG to each 100 mg of gel. For >2% agarose gels, add 6 volumes of Buffer QG. The maximum amount of gel slice per QIAquick column is 400 mg; for gel slices >400 mg use more than one QIAquick column.

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until the gel slice has completely dissolved). To help dissolve gel, mix by vortexing the tube every 2–3 min during the incubation.

**IMPORTANT:** Solubilize agarose completely. For >2% gels, increase incubation time.

4. After the gel slice has dissolved completely, check that the color of the mixture is yellow (similar to Buffer QG without dissolved agarose).

If the color of the mixture is orange or violet, add 10 µl of 3 M sodium acetate, pH 5.0, and mix. The color of the mixture will turn to yellow. The adsorption of DNA to the QIAquick membrane is efficient only at pH ≤7.5. Buffer QG contains a pH indicator which is yellow at pH ≤7.5 and orange or violet at higher pH, allowing easy determination of the optimal pH for DNA binding.

5. Add 1 gel volume of isopropanol to the sample and mix.

For example, if the agarose gel slice is 100 mg, add 100 µl isopropanol. This step increases the yield of DNA fragments <500 bp and >4 kb. For DNA fragments between 500 bp and 4 kb, addition of isopropanol has no effect on yield. Do not centrifuge the sample at this stage.

6. Place a QIAquick spin column in a provided 2 ml collection tube.

7. To bind DNA, apply the sample to the QIAquick column, and centrifuge for 1 min.

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The maximum volume of the column reservoir is 800 µl. For sample volumes of more than 800 µl, simply load and spin again.

**8. Discard flow-through and place QIAquick column back in the same collection tube.**

Collection tubes are reused to reduce plastic waste.

**9. Recommended: Add 0.5 ml of Buffer QG to QIAquick column and centrifuge for 1 min.**

This step will remove all traces of agarose. It is only required when the DNA will subsequently be used for direct sequencing, in vitro transcription, or microinjection.

**10. To wash, add 0.75 ml of Buffer PE to QIAquick column and centrifuge for 1 min.**

**Note:** If the DNA will be used for salt-sensitive applications, such as blunt-end ligation and direct sequencing, let the column stand 2–5 min after addition of Buffer PE, before centrifuging.

**11. Discard the flow-through and centrifuge the QIAquick column for an additional 1 min**

at 17,900 x g (13,000 rpm).

**IMPORTANT:** Residual ethanol from Buffer PE will not be completely removed unless

the flow-through is discarded before this additional centrifugation.

**12. Place QIAquick column into a clean 1.5 ml microcentrifuge tube.**

**13. To elute DNA, add 50 µl of Buffer EB (10 mM Tris·Cl, pH 8.5) or water (pH 7.0–8.5) to the center of the QIAquick membrane and centrifuge the column for 1 min. Alternatively, for increased DNA concentration, add 30 µl elution**

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buffer to the center of the QIAquick membrane, let the column stand for 1 min, and then centrifuge for 1 min.

**IMPORTANT:** Ensure that the elution buffer is dispensed directly onto the QIAquick membrane for complete elution of bound DNA. The average eluate volume is 48 µl from 50 µl elution buffer volume, and 28 µl from 30 µl. Elution efficiency is dependent on pH. The maximum elution efficiency is achieved between pH 7.0 and 8.5. When using water, make sure that the pH value is within this range, and store DNA at -20°C as DNA may degrade in the absence of a buffering agent. The purified DNA can also be eluted in TE (10 mM Tris-Cl, 1 mM EDTA, pH 8.0), but the EDTA may inhibit subsequent enzymatic reactions.

**14. If the purified DNA is to be analyzed on a gel add 1 volume of Loading Dye**

to 5 volumes of purified DNA. Mix the solution by pipetting up and down before loading the gel.

Loading dye contains 3 marker dyes (bromophenol blue, xylene cyanol, and orange G) that facilitate estimation of DNA migration distance and optimization of agarose gel run time. Refer to Table 2 (page 15) to identify the dyes according to migration distance and agarose gel percentage and type.

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In PCR Taq polymerase was used. If Taq polymerase is used, the PCR products include at the 3' end a overhanging A nucleotide. This feature is used by the so-called T/A cloning.

The DNA fragments (SRY gene) are cloned into a vector which has

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overhanging 5' T nucleotides. So these vector endings are compatible to the endings of the PCR products, and the PCR products can be ligated with high efficiency into the vector. One of such vectors is the pGEM-T Easy.

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
Promega pGEM T Easy (Promega Cat. #A1380)

Product	Size	Cat.#
pGEM <sup>®</sup> -T Easy Vector System II	20 reactions	A1380

For Laboratory Use. Includes:

- 1.2µg pGEM<sup>®</sup>-T Easy Vector (50ng/µl)
- 12µl Control Insert DNA (4ng/µl)
- 100u T4 DNA Ligase
- 200µl 2X Rapid Ligation Buffer, T4 DNA Ligase
- 1.2ml JM109 Competent Cells, High Efficiency (6 × 200µl)
- 1 Protocol

Storage Conditions: For Cat.# A3610, A1380, store the Competent Cells at -70°C. All other components can be stored at -20°C.

2X Rapid Ligation Buffer, T4 DNA Ligase (provided with pGEM-T Easy II Kit)

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### 2.3.3 Protocol for Ligations Using the pGEM®-T Easy Vector and the 2X Rapid Ligation Buffer /

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1. Briefly centrifuge the pGEM®-T or pGEM®-T Easy Vector and Control Insert DNA tubes to collect contents at the bottom of the tubes.

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2. Set up ligation reactions as described below.

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**Note:** Use 0.5ml tubes known to have low DNA-binding capacity (e.g., VWR Cat.# 20170-310).

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Vortex the 2X Rapid Ligation Buffer vigorously before each use.

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\*Molar ratio of PCR product:vector may require optimization (see Section VI.C).

3. Mix the reactions by pipetting. Incubate the reactions 1 hour at room temperature. Alternatively, if the maximum number of transformants is required, incubate the reactions overnight at 4°C.

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1. Use only Promega T4 DNA Ligase supplied with this system to perform pGEM®-T and pGEM®-T Easy Vector ligations. Other commercial preparations of T4 DNA ligase may contain exonuclease activities that may remove the terminal . Use only Promega T4 DNA Ligase supplied

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with this system to perform pGEM®-T and pGEM®-T Easy Vector ligations. Other commercial preparations of T4 DNA ligase may contain exonuclease activities that may remove the terminal deoxythymidines from the vector.

2. 2X Rapid Ligation Buffer contains ATP, which degrades during temperature fluctuations. Avoid multiple freeze-thaw cycles and exposure to frequent temperature changes by making single-use aliquots of the buffer.

3. It is important to vortex the 2X Rapid Ligation Buffer before each use.

4. Longer incubation times will increase the number of transformants. Generally, incubation overnight at 4°C will produce the maximum number of transformants from the vector.

2. 2X Rapid Ligation Buffer contains ATP, which degrades during temperature fluctuations. Avoid multiple freeze-thaw cycles and exposure to frequent temperature changes by making single-use aliquots of the buffer.

3. It is important to vortex the 2X Rapid Ligation Buffer before each use.

4. Longer incubation times will increase the number of transformants. Generally, incubation overnight at 4°C will produce the maximum number of transformants.

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For the transformation there are used cells which are already made competent or preperate by the following protocol. These competent E.coli JM109 cells are provided with the Promega pGEM-T Easy Kit II

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1. Set up 5ml culture of E.coli in L.B medium , incubate over night at 37°C on a rotary shaker

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#### 2.4.3 Materials (other than JM 109 from Promega pGEM-T Easy Kit II)<sup>1</sup> / المواد (غير JM 109 from Promega pGEM-T Easy Kit II.)

<sup>1</sup> From Promega, pGEM- T and pGEM-T Easy Systems, Technical Manual, Section XI.C.

• LB plates with ampicillin/IPTG/X-Gal

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LB medium (per liter)  
10g Bacto®-tryptone  
5g Bacto®-yeast extract  
5g NaCl  
Adjust pH to 7.0 with NaOH.

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#### LB plates with ampicillin

Add 15g agar to 1 liter of LB medium. Autoclave. Allow the medium to cool to 50°C before adding ampicillin to a final concentration of 100µg/ml. Pour 30-35ml of medium into 85mm Petri dishes. Let the agar harden. Store at 4°C for up to 1 month or at room temperature for up to 1 week.

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#### LB plates with ampicillin/IPTG/X-Gal

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Make the LB plates with ampicillin as above; then supplement with 0.5mM IPTG and 80µg/ml X-Gal and pour the plates. Alternatively, 100µl of 100mM IPTG and 20µl of 50mg/ml X-Gal may be spread over the surface of an LB-ampicillin plate and allowed to absorb for 30 minutes at 37°C prior to use.

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#### IPTG stock solution (0.1M)

1.2g IPTG  
Add water to 50ml final volume.  
Filtersterilize and store at 4°C.

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#### X-Gal (2ml)

100mg 5-bromo-4-chloro-3-indolyl-β-D-galactoside  
Dissolve in 2ml N,N'-dimethylformamide. Cover with aluminum foil and store at -20°C.

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The genotype of JM109 is *recA1, endA1, gyrA96, thi, hsdR17 (rK-,mK+), relA1, supE44, Δ(lac-proAB), [F', traD36, proAB, lacIqZΔM15]* (5).

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1. Prepare two LB/ampicillin/IPTG/X-Gal plates for each ligation reaction, plus two plates for determining transformation efficiency (see Section "Transformation Control"). Equilibrate the plates to room temperature prior to plating (Step 10).

2. Centrifuge the tubes containing the ligation reactions to collect contents at the bottom of the tube. Add 2µl of each ligation reaction to a sterile 1.5ml microcentrifuge tube on ice (see Note 1). Set up another tube on ice with 0.1ng uncut plasmid for determination of the transformation efficiency of the competent cells (see Section "Transformation Control").

3. Remove tube(s) of frozen JM109 High Efficiency Competent Cells from storage and place in an ice bath until just thawed (about 5 minutes). Mix the cells by **gently** flicking the tube.

**Note:** Avoid excessive pipetting, as the competent cells are extremely fragile.

4. **Carefully** transfer 50µl of cells into each tube prepared in Step 2 (100µl cells for determination of transformation efficiency).

5. **Gently** flick the tubes to mix and place them on ice for 20 minutes.

6. Heat-shock the cells for 45–50 seconds in a water bath at exactly 42°C (**Do not shake**).

7. Immediately return the tubes to ice for 2 minutes.

8. Add 950µl room temperature SOC medium to the tubes containing cells transformed with ligation reactions and 900µl to the tube containing cells transformed with uncut plasmid (LB broth may be substituted, but colony number may be lower).<sup>2</sup>

9. Incubate for 1.5 hours at 37°C with shaking (~150rpm).

10. Plate 100µl of each transformation culture onto duplicate LB/ampicillin/IPTG/X-Gal plates. For the transformation control, a 1:10 dilution with SOC medium is

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recommended for plating. If a higher number of colonies is desired, the cells may be pelleted by centrifugation at  $1,000 \times g$  for 10 minutes, resuspended in  $200\mu\text{l}$  of SOC medium, and  $100\mu\text{l}$  plated on each of two plates.

11. Incubate the plates overnight (16–24 hours) at  $37^\circ\text{C}$ . In our experience, if  $100\mu\text{l}$  is plated approximately 100 colonies per plate are routinely seen when using competent cells that are  $1 \times 10^8 \text{cfu}/\mu\text{g}$  DNA. Longer incubations or storage of plates at  $4^\circ\text{C}$  (after  $37^\circ\text{C}$  overnight incubation) may be used to facilitate blue color development. White colonies generally contain inserts; however, inserts may also be present in blue colonies.

Please see Section “Screening Transformants for Inserts: Preselection of transformants” for more information.

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Notes:

1. In our experience, the use of larger ( $17 \times 100\text{mm}$ ) polypropylene tubes (e.g., Falcon Cat.# 2059) has been observed to increase transformation efficiency. Tubes from some manufacturers bind DNA, thereby decreasing the colony number, and should be avoided.

2. Colonies containing  $\beta$ -galactosidase activity may grow poorly relative to cells lacking this activity. After overnight growth, the blue colonies may be smaller than the white colonies, which are approximately one millimeter in diameter.

Promega recommends to transformation control.<sup>3</sup> This won't be done in this course.

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Successful cloning of an insert into the pGEM®-T and pGEM®-T Easy Vectors interrupts the coding sequence of  $\beta$ -galactosidase; recombinant clones can

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
usually be identified by color screening on indicator plates. However, the characteristics of PCR products cloned into the pGEM®-T and pGEM®-T Easy Vectors can significantly affect the ratio of blue:white colonies obtained following transformation of competent cells. Clones that contain PCR products, in most cases, produce white colonies, but blue colonies can result from PCR fragments that are cloned in-frame with the *lacZ* gene. Such fragments are usually a multiple of 3 base pairs long (including the 3'-A overhangs), and do not contain in-frame stop codons. There have been reports of DNA fragments of up to 2kb that have been cloned in-frame and have produced blue colonies.

Even if your PCR product is not a multiple of 3 bases long, the amplification process can introduce mutations (e.g., deletions or point mutations) that may result in blue colonies when competent cells are transformed with the fragment inserted into the pGEM®-T or pGEM®-T Easy Vectors.

The white colonies have to be verified with plasmid miniprep.

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Set up mini cultures of 4 white colonies in 2.5ml LB/Ampicilin. The cultures are incubated over night at 37°C.

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
Figure: Taking E.coli with sterile soft loop /

## 2.6 Isolation of recombinant plasmid from E.coli cells with Qiagen Miniprep kit (5<sup>th</sup> day)<sup>4</sup> / عزل البلازميد /

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The plasmid has to be isolated from the overnight cultures

QIAprep Spin Miniprep Kit (50)	For 50 high-purity plasmid minipreps: 50	27104
	QIAprep Spin Columns, Reagents, Buffers, Collection Tubes (2 ml)	

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**Storage:** QIAprep Miniprep Kits should be stored dry at room temperature (15–25°C). Kits can be stored for up to 12 months without showing any reduction in performance and quality. For longer storage these kits can be kept at 2–8°C.

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If any precipitate forms in the buffers after storage at 2–8°C it should be redissolved by warming the buffers to 37°C before use. After addition of RNase A and optional LyseBlue reagent, Buffer P1 is stable for 6 months when stored at 2–8°C. RNase A stock solution can be stored for two years at room

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The QIAprep miniprep procedure is based on alkaline lysis of bacterial cells followed by adsorption of DNA onto silica in the presence of high salt (1). The unique silica membrane used in QIAprep Miniprep Kits completely replaces glass or silica slurries for plasmid minipreps. The procedure consists of three basic steps:

■□Preparation and clearing of a bacterial lysate

■□Adsorption of DNA onto the QIAprep membrane

■□Washing and elution of plasmid DNA

All steps are performed without the use of phenol, chloroform, CsCl, ethidium bromide, and without alcohol precipitation.

#### Preparation and clearing of bacterial lysate

The QIAprep miniprep procedure uses the modified alkaline lysis method of Birnboim and Doly (2). Bacteria are lysed under alkaline conditions, and the lysate is subsequently neutralized and adjusted to high-salt binding conditions in one step. After lysate clearing, the sample is ready for purification on the QIAprep silica membrane.

#### DNA adsorption to the QIAprep membrane

QIAprep columns, strips, and plates use a silica membrane for selective adsorption of plasmid DNA in high-salt buffer and elution in low-salt buffer. The optimized buffers in the lysis procedure, combined with the unique silica membrane, ensure that only DNA will be adsorbed, while RNA, cellular proteins, and metabolites are not retained on the membrane but are found in the flow-through.

#### Washing and elution of plasmid DNA

Endonucleases are efficiently removed by a brief wash step with Buffer PB. This step is essential when working with endA+ strains such as the JM series, HB101 and its derivatives, or any wild-type strain, to ensure that plasmid DNA is not degraded.

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The Buffer PB wash step is also necessary when purifying low-copy plasmids, where large culture volumes are used. Salts are efficiently removed by a brief wash step with Buffer PE. High-quality plasmid DNA is then eluted from the QIAprep column with 50–100 µl of Buffer EB or water. The purified DNA is ready for immediate use in a range of applications — no need to precipitate, concentrate, or desalt.

**Note:** Elution efficiency is dependent on pH. The maximum elution efficiency is achieved between pH 7.0 and 8.5. When using water for elution, make sure that the pH value is within this range. Store DNA at -20°C when eluted with water since DNA may degrade in the absence of a buffering agent.

#### DNA yield

Plasmid yield with the QIAprep miniprep system varies depending on plasmid copy number per cell (see page 39), the individual insert in a plasmid, factors that affect growth of the bacterial culture (see pages 39–42), the elution volume (Figure 1), and the elution incubation time (Figure 2). A 1.5 ml overnight culture can yield from 5 to 15 µg of plasmid DNA (Qiagen Miniprep Handbook, Table 1, page 14). To obtain the optimum combination of DNA quality, yield, and concentration, we recommend using Luria-Bertani (LB) medium for growth of cultures (for composition see Qiagen Miniprep Handbook, page 41), eluting plasmid DNA in a volume of 50 µl, and performing a short incubation after addition of the elution buffer.

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
In the following text is described the Plasmid DNA Purification using the QIAprep Spin Miniprep Kit and a Microcentrifuge. This protocol is designed for purification of up to 20 µg of high-copy plasmid DNA from 1–5 ml overnight cultures of E. coli in LB (Luria-Bertani) medium. For purification of low-copy plasmids and cosmids, large plasmids (>10 kb), and DNA

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بروتوكول آخر في كتيب prepared using other methods, refer to the recommendations on Qiagen Miniprep Handbook, page 44. Please read "Important Notes" of Qiagen Miniprep Handbook, pages 15–21 before starting. Note: All protocol steps should be carried out at room temperatur

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1. Resuspend pelleted bacterial cells in 250 µl Buffer P1 and transfer to a microcentrifuge tube.

Ensure that RNase A has been added to Buffer P1. No cell clumps should be visible after resuspension of the pellet. If LyseBlue reagent has been added to Buffer P1, vigorously shake the buffer bottle to ensure LyseBlue particles are completely dissolved. The bacteria should be resuspended completely by vortexing or pipetting up and down until no cell clumps remain.

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by inverting the tube 4–6 times.

Mix gently by inverting the tube. Do not vortex, as this will result in shearing of genomic DNA. If necessary, continue inverting the tube until the solution becomes viscous and slightly clear. Do not allow the lysis reaction to proceed for more than 5 min. If LyseBlue has been added to Buffer P1 the cell suspension will turn blue after addition of Buffer P2. Mixing should result in a homogeneously colored suspension. If the suspension contains localized colorless regions or if brownish cell clumps are still visible, continue mixing the solution until a homogeneously colored suspension is achieved.

3. Add 350 µl Buffer N3 and mix immediately and thoroughly by inverting the tube 4–6 times.

To avoid localized precipitation, mix the solution thoroughly, immediately after addition of Buffer N3. Large culture volumes (e.g. ≥5 ml) may require inverting up to 10 times. The solution should become cloudy. If LyseBlue reagent has been used, the suspension should be mixed until all trace of blue has gone and the suspension is colorless. A homogeneous colorless suspension indicates that the SDS has been effectively precipitated.

4. Centrifuge for 10 min at 13,000 rpm (~17,900 × g) in a table-top microcentrifuge. A compact white pellet will form.

5. Apply the supernatants from step 4 to the QIAprep spin column by decanting or pipetting.

6. Centrifuge for 30–60 s. Discard the flow-through.

7. Recommended: Wash the QIAprep spin column by adding 0.5 ml Buffer PB and centrifuging for 30–60 s. Discard the flow-through.

This step is necessary to remove trace nuclease activity when using endA<sup>+</sup> strains such as the JM series, HB101 and its derivatives, or any wild-type strain, which have high levels of nuclease activity or high carbohydrate content. Host strains such as XL-1 Blue and DH5<sup>TM</sup> do not require this

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9. Discard the flow-through, and centrifuge for an additional 1 min to remove residual wash buffer.

Important: Residual wash buffer will not be completely removed unless the flow-through is discarded before this additional centrifugation. Residual ethanol from Buffer PE may inhibit subsequent enzymatic reactions.

10. Place the QIAprep column in a clean 1.5 ml microcentrifuge tube. To elute DNA, add 50 µl Buffer EB (10 mM Tris·Cl, pH 8.5) or water to the center of each QIAprep spin column, let stand for 1 min, and centrifuge for 1 min.

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With the restriction digest we want to confirm the a proper DNA fragment has been cloned into the plasmid vector.

On both sides of the insertion site of pGEM-T Easy is a EcoRI interface, that means that the insert could be cut completely out of the plasmid with EcoRI digest.

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1. Set up restriction digest and incubate for at least 1 hour at 37°C.

Restriction digest set up:

1µl Plasmid DNA

2µl 10xBuffer

0.2µl EcoRI

16.8 µl H<sub>2</sub>O (ad 20 µl with H<sub>2</sub>O)

20µl

2. After 1 hour of digestion 2 µl

probe buffer is added and the whole probe is put onto a 1% agarose gel. Also add a DNA

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ladder into another column to be able to identify the largeness of analysed DNA molecule. There should be a DNA with 897bp.

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There are two tests for the new influenza virus (causing the «swine influenza»):

- Fast antibody test
- Genetic (PCR) test

The fast antibody test is not reliable, because it has a failing rate of about 50%.

The Genetic (PCR) test is today done with a real time (RT) PCR machine. This machine costs at least about 40.000 USD. Here is no Gel electrophoresis needed.

If no RT PCR machine is available, but only a PCR machine, the test could be done by nested PCR.

The key to successful purification of RNA from cells and tissues is speed . Cellular Rnases should be inactivated as quickly as possible at the very first stage in the extraction process . Once the endogenous Rnases have been destroyed , the immediate threat to the integrity of the RNA is greatly reduced , and purification can proceed at a more graceful pace.

Because of the urgency ,many methods for the isolation of the intact RNA from cells use strong denaturants such as guanidinium hydrochloride or guanidinium thiocyanate to disrupt cells , solubilize their components, and denature endogenous Rnases simultaneously. The use of guanidinium isocyanate in RNA extraction, first mentioned briefly by Ullrich et al.1977 , was documented in papers published by Han et al. 1987 and Chirgwin et al . 1979. The Han method is laborious as it involves solubilization of RNA pellets in progressively smaller volumes of 5 M guanidinium thiocyanate. In the Chirgwin method, cultured cells or tissues are homogenized in 4 M guanidinium isothiocyanate, and the

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lysate is layered onto a dense cushion of CsCl . Because the buoyant density of RNA in CsCl (1.8 g/ml ) is much greater than that of other cellular components, rRNAs and mRNAs migrate to the bottom of the tube during ultracentrifugation (Glisin et al. 1974 ). As long as the step gradients are not overloaded, proteins remain in the guanidium lysate while DAN floats on the CsCl cushion. Because the Chirgwin method yields RNA of very high quality and purity and is not labor/intensive, it became the standard technique during the early 1980s for isolation of undergraded high molecular weight RNA . However, the method has one weakness; IT is unsuitable for simultaneous processing of many samples. For this purpose, it has been almost completely displaced by the single step technique of Chomczynski and Sacchi ( 1987), in which the guanidinium thiocyanate homogenate is extracted of with phenol;chloroform at reduced PH. Elimination of the ultracentrifugation step allows many samples to be processe simultaneously and speedily at modest cost and without sacrifice in yield or quality of RNA. For many investigatorse, the single step technique described in protocol 1 remains the method of choice to isolate RNA from cultured cells and most animal tissues.

There are two circumstances in which single-step procedure is not recommended. First , the procedure does not extract RNA efficiently from adipose tissues that are rich in triglycerides. RNA is best prepared from these fatty sources by a modification of the method of Chirgwin et al. , described by Tavangar et al . Second, RNA prepared by guanidine lysis is sometimes contaminated to a

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significant extent by cellular polysaccharides and proteoglycans . These contaminant are reported to prevent solubilization of RNA after precipitation with alcohols , to inhibit Reverse-transcriptase-polymerase chain reaction (RT-PCRs), and to bind to membranes during RNA blotting . If contamination by proteoglycans appears to be a problem , include an organic extraction step and change the condition used to precipitate the RNA .

The yield of total RNA depends on the tissue or cell source, but it is generally in the range of 4-7µg/mg of starting tissue or 5-10 µg/10<sup>6</sup> cells.

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This protocol describes a single-step for the purification of RNA. Cells are homogenized in guanidinium thiocyanate and the RNA is purified from the lysate by extraction with phenol-chloroform at reduced PH. This method allows many samples to be processed simultaneously and quickly. The yield of total RNA depends on the tissue or cell source and is generally in the range of 4-7µg per ml starting tissue or 5-10µg per 10<sup>6</sup> cells. All reagents used in this protocol must be prepared with diethyl pyrocarbonate (DEPC)-treated H<sub>2</sub>O.

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#### Procedure:

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1. Prepare cells or tissue samples for isolation of RNA as appropriate for the material under study. Consult the table below for the amounts of the solution D (4 M guanidinium thiocyanate , 25 mM sodium citrate-2H<sub>2</sub>O, 0.5% (wt/vol) sodium lauryl sarcosinate, 0.1 M 2-mercaptoethanol) required for different types of samples.

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#### Amount of tissue or cells / Amount of solution D

100mg of tissue	3 ml
75ml(1-75)flask of cells	3 ml
60mm plate of cells	1ml
90mm plate of cells	2 ml

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#### For tissues:

A-Isolate the desired tissues by dissection and place them immediately in liquid nitrogen.

B-Transfer 100 mg of the frozen tissue to a mortar containing liquid nitrogen and pulverize the tissue using a pestle. Keep the tissue frozen by the addition of liquid nitrogen.

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polpropylene snap-cap tube containing 3 ml of solution D and homogenized at 15-25 °C for 15-30 s with a polytron homogenizer(kinematica).

For mammalian cells grown in suspension:

A-Harvest the cells by centrifugation at 200-1900g at 15-25 °C for 5-10 min. Resuspend the cells in 1-2 ml of sterile ice-cold PBS (137mM NaCl, 2.7mM KCl, 10mM Na<sub>2</sub>HPO<sub>4</sub>, 2mM KH<sub>2</sub>PO<sub>4</sub>).

B-Harvest the cells again by centrifugation, remove the PBS completely by aspiration and add 2 ml of solution D per 10<sup>+</sup> cells.

C-Homogenize the lysates with a Polytron homogenizer at 15-25 °C for 15-30 s.

For mammalian cells grown in monolayers:

A-remove the medium and rinse the cells once with 5-10 ml of sterile ice-cold PBS.

B-remove PBS and lyse the cells in 2ml of solution D per 90mm culture dish (1ml per 60mm dish). Transfer the lysates to a polypropylene snap-cap tube.

C-Homogenize the lysates with a polytron homogenizer at 15-25 °C for 15-30 s.

2-Transfer the homogenate to a fresh tube and sequentially add 0.1ml of 2 M sodium acetate (PH 4.0), 1 ml of phenol and 0.2 ml of 49:1(wt/vol) chloroform/ isoamyl per millilitre of solution D. After addition of each reagent, cap the tube and mix the contents thoroughly by inversion.

3-Vortex the homogenate vigorously for 10 s. Incubate the tube for 15 min on ice.

4-Centifuge the tube at 10,000 g at 4 °C for 20 min and then transfer the

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upper, aqueous phase containing the extracted RNA to a fresh tube.

5-Add an equal volume of isopropanol to the extracted RNA. Mix the solution well and allow the RNA to precipitate at  $-20^{\circ}\text{C}$  for at least 1 h.

6-Collect the precipitated RNA by centrifugation at  $10.000\text{g}$  at  $^{\circ}\text{C}$  for 30 min.

7-Carefully decant the isopropanol and dissolve the RNA pellet in 0.3 ml of solution D for every 1ml of this solution used in step 1.

8-Transfer the solution to a microcentrifuge tube, vortex it well and precipitate the RNA with 1 volume of isopropanol at  $-20^{\circ}\text{C}$  for 1 h or more.

Collect the precipitated RNA by centrifugation at maximum speed at  $4^{\circ}\text{C}$  for 10 min in a microcentrifuge. Wash the pellet twice with 75% ethanol, centrifuge again and remove any remaining ethanol with a disposable pipette tip. Allow the pellet to air dry for a few minutes before dissolving it in 50-100  $\mu\text{l}$  of DEPC-treated  $\text{H}_2\text{O}$  and storing at  $-70^{\circ}\text{C}$ .

9-Estimate the concentration of the RNA by measuring the absorbance at 260 nm of an aliquot of the final preparation.

Purified RNA is not immune to degradation by RNase after resuspension in the 0.5% SDS solution. Some investigators therefore prefer to dissolve the pellet of RNA in 50-100  $\mu\text{l}$  of stabilized formamide and store the solution at  $-20^{\circ}\text{C}$ . RNA can be recovered from formamide by precipitation with 4 volumes of ethanol. For further details, please see the panel on STORAGE OF RNA.

SDS should be removed by chloroform extraction and ethanol precipitation before enzymatic treatment of the

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RNA(e.g, primer extension,reverse transcriptase, and in vitro translation). The redissolved RNA can then be used for mRNA purification by oligo (dT)-cellulose chromatography or analysed by standard techniques such as blot hybridization or mapping.

RNA prepared from tissues is generally not contaminated to a significant extent with DNA. However, RNA prepared from cell lines undergoing spontaneous or induced apoptosis is often contaminated with fragments of degraded genomic DNA, RNA prepared from transfected cells is almost always contaminated by fragments of the RNA used for transfection. Some investigators therefore treat the final RNA preparation with RNase-free DNase. Alternatively, fragments of DNA may be removed by preparing poly(A)<sup>+</sup> RNA by oligo(dT)chromatography.

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#### 4.1 STORAGE OF RNA / تخزين الـ RNA

After precipitation with ethanol, store the RNAs as follows:

- Dissolve the precipitate in deionized and store at -20 °C. Formamide provides a chemically stable environment that also protects RNA against degradations by RNases. Purified, salt-free RNA dissolves quickly in formamide up to a concentration of 4 mg/ml. At such concentrations, samples of the RNA can be analysed directly by gel electrophoresis, RT-PCR, or RNase protection, saving time and avoiding potential degradation. If necessary, RNA can be recovered from formamide by precipitation with 4 volumes of ethanol as described by Chomczynski or by diluting the formamide fourfold with 0.2 M NaCl and then adding the conventional 2 volumes of ethanol.
- Dissolve the precipitate in an aqueous buffer and store at -80 °C. Buffers commonly used for this purpose include

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SDS (0.1-0.5%) in TE (PH 7.6) or in DEPC-treated H<sub>2</sub>O containing 0.1 mM EDTA (PH 7.5). The SDS should be removed by chloroform extraction and ethanol precipitation before enzymatic treatment of the RNA (e.g, primer extension, reverse transcription, and in vitro translation).

- Store the precipitate of RNA as a suspension at -20 °C in ethanol. Samples of the RNA can be removed, as needed, with an automatic pipetting device. However, because precipitates of RNA are lumpy and sticky, and partly because of losses onto the surfaces of disposable pipette tips, the recovery of RNA is inconsistent.

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[http://nar.oxfordjournals.org/cgi/pdf\\_extract/19/14/4011](http://nar.oxfordjournals.org/cgi/pdf_extract/19/14/4011); *A Rapid membrane-based viral RNA isolation method for the polymerase chain reaction*, Kevin R. Porter et.al., Naval Medical Research Institute, Bethesda, MD and the Walter Reed Army institute for Research, Washington, DC, USA

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Methods used to prepare RNA for use in PCR generally require the use of hazardous chemicals such as guanidine isothiocyanate, chloroform and phenol.

Here we describe a rapid membrane-based method for the isolation of RNA from viruses present in serum. In this system whole virus particles present in serum are immobilized on a membrane by positive pressure filtration. Digestion of the membrane-bound particle with a detergent and proteinase K releases nucleic acid from its capsid and destroys ribonucleases present in the sample.

Strips of Immobilon-P membrane (polyvinylidene difluoride (PVDF), Millipore Corporation, Bedford, MA), measuring 5/32 inches by 1 1/2 inches, were placed on three (3) thicknesses of Whatman 3MM filter paper. One hundred microliters of human serum spiked with varying amounts of West Nile virus (Family Flaviviridae, RNA genome), strain 956B, were spotted near the end of the membrane in two 50 µl spots. A 1 ml tuberculin syringe was used to apply the samples by positive pressure (1, 2). The membrane strip spots were immediately immersed into 0.5 ml microfuge tubes containing 81.75 µl of digestion buffer (79.25 µl distilled water, 0.5 µl Tween-20 and 2 µl proteinase K (20 mg/ml)) and the tubes closed with a portion of the membrane extending outside. The tubes were incubated at 50°C for one hour and subsequently heated to 95°C for 10 minutes to inactivate the proteinase K. After removal of the membranes, 1 µl (0.33 µg) each of sense and antisense strand primer was added. The specimens were heated to 68°C and cooled on ice

to allow primer annealing. Sequentially added were 5  $\mu$ l 2 mM stock of deoxyribonucleotide triphosphates, 10  $\mu$ l 10 $\times$ PCR buffer, 0.5  $\mu$ l RNasin (Sigma Chemical, St Louis, MO), and 0.25  $\mu$ l (200 units/ $\mu$ l) of murine maloney leukemia virus reverse transcriptase. Following a one hour incubation at 42°C, brief denaturation at 94°C and a second 42°C incubation, 0.5  $\mu$ l of Taq DNA polymerase (1.25 units) was added to give a final volume of 100  $\mu$ l. Samples were then amplified using a PCR protocol for West Nile virus (3).

Ten microliters of PCR product were Southern transferred to Nytran membrane. The membrane was hybridized to a 151-base pair digoxigenin-labeled probe (Genius Kit, Boehringer Mannheim Biochemicals, Indianapolis, IN). Hybridized probe was detected by chemiluminescence using Lumi-Phos (Boehringer Mannheim Corporation, Indianapolis, IN).

A positive signal was seen in lanes A1 through A5 (Figure 1). No signal was seen with human RNA or dengue virus RNA, another flavivirus. Virus-free serum also failed to show any bands.

The data presented show this membrane-based procedure to be rapid, sensitive and less cumbersome than conventional viral RNA isolation techniques. The hydrophobic Immobilon-P membrane has a very high protein binding capacity. Its nucleic acid binding efficiency is extremely low. This allows any nucleic acid that is liberated to stay in solution while much of the protein can be adsorbed back onto the membrane. This pressure blotting membrane technique has been used previously for the detection of antibodies to other togaviruses (1) and meningococcal polysaccharide (2).

This membrane-based method was able to detect as few as 6.5 plaque forming units of RNA virus in 100  $\mu$ l of serum. This degree of sensitivity should be adequate for the detection of RNA viruses that produce low-level viremias.

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PCR is a powerful method to amplify specific sequences of DNA from a large complex mixture of DNA. For example, you can design PCR primers to amplify a single locus from an entire genome. From a single template molecule, you can

produce over 1 billion copies of the PCR product very quickly. However, the capacity to amplify over one billion fold also increases the possibility of amplifying the wrong DNA sequence over one billion times. The specificity of PCR is determined by the specificity of the PCR primers. For example, if your primers bind to more than one locus (e.g. paralog or common domain), then more than one segment of DNA will be amplified. To control for these possibilities, investigators often employ nested primers to ensure specificity.

Nested PCR means that two pairs of PCR primers were used for a single locus (figure 1). The first pair amplified the locus as seen in any PCR experiment. The second pair of primers (**nested primers**) bind within the first PCR product (figure 4) and produce a second PCR product that will be shorter than the first one (figure 5). The logic behind this strategy is that if the wrong locus were amplified by mistake, the probability is very low that it would also be amplified a second time by a second pair of primers.

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**Figure 1.** Nested PCR strategy. Segment of DNA with dots representing nondiscript DNA sequence of unspecified length. The double lines represent a large distance between the portion of DNA illustrated in this figure. The portions

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of DNA shown with four bases in a row represent PCR primer binding sites, though real primers would be longer.


**Figure 2.** The first pair of PCR primers (blue with arrows) bind to the outer pair of primer binding sites and amplify all the DNA in between these two sites.

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**Figure 3.** PCR product after the first round of amplification. Notice that the bases outside the PCR primer pair are not present in the product.

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**Figure 4.** Second pair of **nested primers** (red with arrows) bind to the first PCR product. The binding sites for the second pair of primers are a few bases "internal" to the first primer binding sites.

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**Figure 5.** Final PCR product after second round of PCR. The length of the product is defined by the location of the internal primer binding sites.

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When a complete genome sequence is known, it is easier to be sure you will not amplify the wrong locus but since very few of the world's genomes have been sequenced completely, nested primers will continue to be an important control for many experiments.

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agaaatgctggatctggtattatcatttcagatacaccagtcacgattgcaatacaactgtcagacaccaag

ggtgctataaacaccagcctcccatttcagaatatacatccgatcacaattggaaaatgt

Remarque: before  
starting put a  
FILTOSTAT FS  
CODE

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Protocol: Purification of Viral RNA (Spin Protocol)

This protocol is for purification of viral RNA from 140 µl plasma, serum, urine, cellculture

#### Important points before starting

■ Read "Important Notes" (pages-98)

■ All centrifugation steps are carried out at room temperature (15–25°C).

#### Things to do before starting

■ Equilibrate samples to room temperature (15–25°C).

■ Equilibrate Buffer AVE to room temperature for elution in step 11.

■ Check that Buffer AW1 and Buffer AW2 have been prepared according to the instructions at the end of this chapter

■ Add carrier RNA reconstituted in Buffer AVE to Buffer AVL

**Important:** if not have QIAamp Viral RNA mini , we can use peqGOLD Viral RNA Kit protocols

See page 33 or visit [www.peqlab.de](http://www.peqlab.de)

#### Procedure

1. Pipet 560 µl of prepared Buffer AVL containing carrier RNA into a 1.5 ml microcentrifuge tube.

If the sample volume is larger than 140 µl, increase the amount of

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Buffer AVL-carrier RNA proportionally (e.g., a 280 µl sample will require 1120 µl Buffer AVL-carrier RNA) and use a larger tube.

\* Fully automatable on the QIAcube. See [www.qiagen.com/MyQIAcube](http://www.qiagen.com/MyQIAcube) for protocols.

**2. Add 140 µl plasma, serum, urine, cell-culture supernatant, or cell-free body fluid to the Buffer AVL-carrier RNA in the microcentrifuge tube. Mix by pulse-vortexing for 15 s.**

To ensure efficient lysis, it is essential that the sample is mixed thoroughly with Buffer AVL to yield a homogeneous solution. Frozen samples that have only been thawed once can also be used.

**3. Incubate at room temperature (15–25°C) for 10 min.**

Viral particle lysis is complete after lysis for 10 min at room temperature. Longer incubation times have no effect on the yield or quality of the purified RNA.

Potentially infectious agents and RNases are inactivated in Buffer AVL.

**4. Briefly centrifuge the tube to remove drops from the inside of the lid.**

**5. Add 560 µl of ethanol (96–100%) to the sample, and mix by pulse-vortexing for 15 s. After mixing, briefly centrifuge the tube to remove drops from inside the lid.**

Only ethanol should be used since other alcohols may result in reduced RNA yield and purity. Do not use denatured alcohol, which contains other substances such as methanol or methylethylketone. If the sample volume is greater than 140 µl, increase the amount of ethanol proportionally (e.g., a 280 µl sample will require 1120 µl of ethanol). In order to ensure efficient binding, it is essential that the sample is mixed thoroughly with the ethanol to yield

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a homogeneous solution.

**6. Carefully apply 630 µl of the solution from step 5 to the QIAamp Mini column (in a 2 ml collection tube) without wetting the rim. Close the cap, and centrifuge at 6000 x g (8000 rpm) for 1 min. Place the QIAamp Mini column into a clean 2 ml collection tube, and discard the tube containing the filtrate.**

Close each spin column in order to avoid cross-contamination during centrifugation. Centrifugation is performed at 6000 x g (8000 rpm) in order to limit microcentrifuge noise. Centrifugation at full speed will not affect the yield or purity of the viral RNA. If the solution has not completely passed through the membrane, centrifuge again at a higher speed until all of the solution has passed through.

**7. Carefully open the QIAamp Mini column, and repeat step 6.**

If the sample volume was greater than 140 µl, repeat this step until all of the lysate has been loaded onto the spin column.

**8. Carefully open the QIAamp Mini column, and add 500 µl of Buffer AW1. Close the cap, and centrifuge at 6000 x g (8000 rpm) for 1 min. Place the QIAamp Mini column in a clean 2 ml collection tube (provided), and discard the tube containing the filtrate.**

It is not necessary to increase the volume of Buffer AW1 even if the original sample volume was larger than 140 µl.

**9. Carefully open the QIAamp Mini column, and add 500 µl of Buffer AW2. Close the cap and centrifuge at full speed (20,000 x g; 14,000 rpm) for 3 min. Continue directly with step 11, or to eliminate any chance of possible Buffer AW2 carryover, perform step 10, and then continue with**

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**step 11.**

**Note:** Residual Buffer AW2 in the eluate may cause problems in downstream applications. Some centrifuge rotors may vibrate upon deceleration, resulting in flow-through, containing Buffer AW2, contacting the QIAamp Mini column. Removing the QIAamp Mini column and collection tube from the rotor may also cause flow-through to come into contact with the QIAamp Mini column. In these cases, the optional step 10 should be performed.

**10. Recommended:** Place the QIAamp Mini column in a new 2 ml collection tube (not provided), and discard the old collection tube with the filtrate. Centrifuge at full speed for 1 min.

**11. Place the QIAamp Mini column in a clean 1.5 ml microcentrifuge tube (not provided). Discard the old collection tube containing the filtrate. Carefully open the QIAamp Mini column and add 60 µl of Buffer AVE equilibrated to room temperature. Close the cap, and incubate at room temperature for 1 min. Centrifuge at 6000 x g (8000 rpm) for 1 min.**

A single elution with 60 µl of Buffer AVE is sufficient to elute at least 90% of the viral RNA from the QIAamp Mini column. Performing a double elution using 2 x 40 µl of Buffer AVE will increase yield by up to 10%. Elution with volumes of less than 30 µl will lead to reduced yields and will not increase the final concentration of RNA in the eluate.

Viral RNA is stable for up to one year when stored at -20°C or -70°C.

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First-strand cDNA Synthesis Using M-MLV Reverse Transcriptase /

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A 20µl reaction volume can be used for 1ng-5µg of total RNA or 1-500ng of mRNA

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2/-Heat mixture to 65oC for 5 min and quick chill on ice.

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6- Inactivate the reaction by heating at 70°C for 15 min.

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8- Immediately put the tube into ice for 1 min

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11. Add the above to a PCR reaction tube for a final reaction volume of 50 µl :

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

7- take a image with digital camera (with a red filter)

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Result: The person is infected with H1N1 virus because there are bands seen /

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If preparing RNA for the first time please read "Handling RNA" (page 34). All steps of the QIAamp Viral RNA Mini protocols should be performed quickly and at room temperature.

After collection and centrifugation, plasma (untreated or treated with anticoagulants other than heparin) or serum can be stored at 2-8°C for up to 6 hours. For long-term storage, freezing at -20°C to -80°C in aliquots is recommended. Frozen plasma or serum samples must not be thawed more than once. Repeated freezing and thawing leads to denaturation and precipitation of proteins, causing reduced viral titers and

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subsequently reduced yields of the isolated viral RNA. In addition, cryoprecipitates formed by freeze-thawing will cause clogging of the QIAamp membrane. If cryoprecipitates are visible, they can be pelleted by briefly centrifuging at 6800 x g for 3 minutes. The cleared supernatant should be removed, without disturbing the pellet, and processed immediately. This step will not reduce viral titers.

The QIAamp Viral RNA Mini procedure is not designed to separate RNA from DNA.

To avoid cellular DNA contamination follow the guidelines in "Cellular DNA". The QIAamp Viral RNA Mini procedure isolates all RNA molecules larger than 200 nucleotides. Smaller RNA molecules will not bind quantitatively under the conditions used.

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Ribonucleases (RNases) are very stable and active enzymes that generally do not require cofactors to function. Since RNases are difficult to inactivate and only minute amounts are sufficient to destroy RNA, do not use any plasticware or glassware without first eliminating possible RNase contamination. Great care should be taken to avoid inadvertently introducing RNases into the RNA sample during or after the purification procedure. In order to create and maintain an RNase-free environment, the following precautions must be taken during pretreatment and use of disposable and non-disposable vessels and solutions while working with RNA.

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concentrate. Before using for the first time, add the appropriate amount of ethanol (96-100%) as indicated on the bottle in Buffer AW1 is stable for 1 year when stored closed at room temperature, but only until the kit expiration date.

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Buffer AW2 is supplied as a concentrate. Before using for the first time, add the appropriate amount of ethanol (96-100%) to Buffer AW2 concentrate as indicated on the bottle.

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Buffer AW2 is stable for 1 year when stored closed at room temperature, but only until the kit expiration date.

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Materials required, but not supplied:

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1. Lysis

Add 450 µl RNA Lysis Buffer T to an Extraction Tube. Pipet 150 µl plasma, cell free body fluid, cell culture supernatant or urine into the tube, mix thoroughly by vortexing for 10 seconds and incubate for 15 minutes at room temperature.

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2. Load and bind

Add 600 µl RNA Binding Solution and mix thoroughly by pipetting until a homogeneous solution is formed. Apply 600 µl of the sample to a PerfectBind RNA Column assembled in a 2 ml Collection Tube. Centrifuge for 1 minute at 10.000 x g. Discard Collection Tube and

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assemble Perfect Bind RNA Column to a new Collection Tube. Apply residual 600 µl of the sample and centrifuge for 1 minute at 10.000 x g. Discard Collection Tube and assemble PerfectBind RNA Column to a new Collection Tube.

Note: If the sample should not have run by the filter completely, centrifugation time might be extended!

### 3. Wash I

Add 500 µl RNA Wash Buffer I to the column and centrifuge the PerfectBind RNA

Column/Collection Tube assembly for 1 minute at 10.000 x g. Discard flow-throw liquid and re-use the Collection Tube.

4. DNase I Digestion (optional) Since PerfectBind RNA resin and spin-column technology actually removes most of DNA without the DNase treatment, it is not necessary to do DNase digestion for most downstream applications. However, certain sensitive RNA applications might require further DNA removal. Following steps provide on-membrane DNase I digestion (Order No. 12-1091).

a. For each PerfectBind RNA Column, prepare this DNase I digestion reaction mix:

DNase I Digestion Buffer 73.5 µl  
RNase-free DNase I (20 Kunitz units/µl) 1.5 µl

Total volume 75 µl

Note:

1. DNase I is very sensitive for physical denaturation, so do not vortex this DNase I mixture!

Mix gently by inverting the tube. Prepare the fresh DNase I digestion mixture directly before RNA isolation.

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supplied with RNase-free DNase set. Standard DNase buffers are not compatible with on-membrane DNase digestion!

b. Pipet 75 µl of the DNase I digestion reaction mix directly onto the surface of PerfectBind RNA resin in each column. Make sure to pipet the DNase I digestion mixture directly onto the membrane. DNase I digestion will not be complete if some of the mix stick to the wall or the O-ring of the PerfectBind RNA Column.

c. Incubate at room temperature (25 - 30 °C) for 15 minutes.

d. Place a PerfectBind RNA Column into a new 2 ml Collection Tube and add 400 µl RNA Wash Buffer I. Place the column at benchtop for 5 minutes. Centrifuge at 10.000 x g for 5 minutes and discard flow-through. Re-use Collection Tube in the next step.

Continue with step 5.

#### 5. Wash II

Add 650 µl completed RNA Wash Buffer II to the column and centrifuge the PerfectBind RNA Column/Collection Tube assembly for 1 minute at 10.000 x g. Discard the flowthrough liquid. Repeat this wash step using the same Collection Tube and discard the flow-through liquid.

6. Dry (Important, do not skip this step!)

Place the PerfectBind RNA Column in the Collection Tube and centrifuge for 2 minutes at 10.000 x g to dry the column matrix.

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80 µl RNase-free Water directly to the binding matrix in the PerfectBind RNA Column, incubate 2 minutes and centrifuge for 1 minute at 6.000 x g to elute RNA.

A second elution may be necessary if the expected yield of RNA is > 50 µg. Alternatively, RNA may be eluted with a higher volume of water. While additional elution increase total RNA yield, the concentration will be lowered since more than 80 % of RNA is recovered with the first elution.

Pre-heating RNase-free Water to 70 °C before adding to the spin column and incubating the spin column for 5 minutes at room temperature before centrifugation may increase yield. Instead of eluting twice, the RNA can directly be eluted in a bigger volume of RNase-free Water.

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1 day training course for researchers and technical staff

Authors:

Noha Abdulwahab, Ghina al-Eter, Omran Zaki, Layal Chbib, and Mirna Khoder

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Chickens are susceptible to many infectious diseases. One of the most important of these is the viral disease known as influenza which is caused for examples by the actual H1N1 virus. For this reasons viruses can be propagated will in chicken eggs for vaccine production purpose.

Eggs for our work can be purchase from normal chicken farms.the eggs must be 9-10 days old when purchased.

Then the virus probe is inoculated into the eggs where the virus is propagated while the eggs are incubated. The work must be done under very clean circumstances and atmosphere to prevent contamination.

Testing the actual state of the eggs is done by candling.

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Laboratory staff should be familiar with and have practiced the following skills prior to the commencement of H1N1 disease vaccine production - this manual does not contain further details about these skills:

- Aseptic technique.
- Sterilization by autoclaving and hot air of glassware and discarded materials.

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An order can be placed for the delivery of the eggs. It is useful if the person responsible for placing the orders and receiving the eggs keeps records. The following information should be recorded in a notebook set aside for this purpose.

- Date when the eggs are ordered and the name of the person who received the order.
- Number and age of the eggs ordered.
- Date and number of the eggs received.
- Colour and appearance of the eggs received.
- Number of eggs damaged during

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- transport.
- Date and number of eggs placed in incubator.
- Number of viable eggs after candling prior to inoculation.

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Appropriate chemical disinfectants must be used for cleaning equipment, materials and work surfaces. All laboratory wastes must be assigned to a category and placed in clearly labeled bins from where the waste will be disposed of appropriately.

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#### Alcohol

A 70% volume/volume (v/v) solution of alcohol diluted with water is useful for wiping down benches and disinfecting the outside of eggs before inoculation and harvesting of allantoic fluid. The addition of 2 percent iodine will increase the effectiveness of this solution.

Note that 70 percent solution of alcohol is flammable!

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#### Chlorine

There are several chlorine compounds that are used as disinfectants. Sodium hypochlorite (NaOCl) or household bleach is readily available and cheap.

Soaking overnight in a 2 percent solution of chlorine is useful for disinfecting plastic materials. Note that commercial bleach contains 12 to 14 percent hypochlorite when manufactured but this concentration deteriorates with time.

Note that chlorine damages fabric and corrodes many metals!

Always read the instructions before using disinfectants and cleaning reagents!

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Many vaccine production centres will already have large commercial incubators installed. Smaller incubators are available

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**Note:** In some eggs the air sac will have not developed on the blunt end but half way down the egg. These eggs are not suitable for vaccine production.

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The most convenient method of propagating of influenza virus in the laboratory is by the inoculation of the allantoic cavity of embryonated eggs. All strains of virus will grow in the cell of the allantoic cavity. The virus enters these cells where it multiplies. As the cells are disrupted the virus is shed into the allantoic fluid. Virulent strains of the virus will invade cells beyond the lining of the allantoic cavity and kill the embryo. The time taken for this to occur is the basis of the "Mean Death Time Assays", which indicates the level of virulence. The avirulent strain of influenza virus will not kill embryos inoculated into the allantoic cavity.

Inoculation of the allantoic cavity of embryonated eggs is a technique used in the following procedures:

1. influenza virus vaccine production
2. Establishing the infectivity titre of a suspension of virus.
3. Isolation of virus from field specimens for laboratory Diagnosis

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### Materials

Eggs 9-day old or 10-day old embryonated eggs. Candle the eggs and mark the inoculation sites . Eggs should be placed in an egg rack with the inoculation site uppermost.

- Egg shell punch or forceps.
- Cotton.
- A 70 % alcohol solution in water.
- Syringe 1 mL.
- Needles preferably 25 gauge, 16 mm.
- sticky tape or melted wax to seal the inoculation site.
- Inoculum. This must be free of microbial contamination.
- Discard tray.

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### Method

1. Use cotton wool and 70 percent alcohol to swab the end of the eggs to be inoculated. Allow the alcohol to evaporate.
2. Swab the eggshell punch with 70 percent alcohol solution. Place used cotton wool in discard tray.

3. Pierce a hole in the end of the egg at the marked inoculation site.
4. Attach needle to 1 mL syringe.
5. Draw inoculum into 1 mL syringe.
6. Keeping the needle and syringe vertical, place the needle through the hole in the eggshell. The needle will need to penetrate approximately 16mm into the egg to reach the allantoic cavity.
7. Inject 0.1 mL of inoculum into the egg.
8. Withdraw the needle from the egg.
9. Seal the hole in the shell with stationery tape or melted wax.
10. Discard the used needles and syringes.
11. Place the inoculated eggs into a second incubator. Check the temperature and humidity of incubator.


Figure 2: Inoculation of the allantoic cavity

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The following method describes harvesting a small sample of allantoic fluid for testing for the presence of

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haemagglutinin using the rapid or micro tests.

#### Materials

- Forceps or a small pair of scissors
- Absolute alcohol for flaming forceps
- Cotton.
- 70 percent alcohol solution
- Discard tray
- 50 µL micropipette and tips, a wire loop or sterile Pasteur pipettes

#### Method

1. Chill eggs at 4°C for at least two hours to kill the embryo and to reduce the contamination of the allantoic fluid with blood during harvesting.
2. Remove sticky tape (if used to seal the eggs) and swab each egg with cotton wool soaked with 70 percent alcohol to disinfect and remove condensation from the shells.
3. Dip the forceps or scissors in absolute alcohol and flame to sterilize. Remove the eggshell above the air space.
4. Discard embryos that are visibly contaminated.
5. Remove a sample of allantoic fluid from each egg. Use a micropipette and sterile tip, sterile glass pipette or a flamed loop and dispense the sample according to method being used for the test

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This is the result of the haemagglutinin part of the haemagglutinin/ neuraminidase viral protein binding to receptors on the membrane of red blood cells. The linking together of the red blood cells by the viral particles results in clumping. This clumping is known as haemagglutination. Haemagglutination is visible macroscopically and is the basis of haemagglutination tests to detect the presence of viral particles. The test does not discriminate between viral particles that are infectious and particles that are degraded and no longer able to infect cells. Both can cause the agglutination of red blood cells.

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Every time a haemagglutination test is carried out, it is necessary to test the settling pattern of the suspension of red blood cells. This involves mixing diluent with red blood cells and allowing the cells to settle.

1. Dispense diluent.
2. Add red blood cells and mix by gently shaking.
3. Allow the red blood cells to settle and observe the pattern.
4. Observe if the cells have a normal settling pattern and there is no auto-agglutination. This will be a distinct button of cells in the micro test and an even suspension with no signs of clumping in the rapid test.

**Note:** The diluent used for haemagglutination tests in this manual is PBS. There should be no signs of haemolysis in the red blood cell suspension. If there are signs of haemolysis, a fresh suspension must be prepared. There should not be any sign of auto-agglutination in the red blood cell control. If an agglutination pattern is observed, discard the suspension of red blood cells. Prepare a fresh suspension and test again.

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Tissue culture is often a generic term that refers to both organ culture and cell culture and the terms are often used interchangeably. Cell cultures are derived from either primary tissue explants or cell suspensions. Primary cell cultures typically will have a finite life span in culture whereas continuous cell lines are, by definition, abnormal and are often transformed cell lines.

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When cells are taken freshly from animal tissue and placed in culture, the cultures consist of a wide variety of cell types, most of which are capable of very limited growth in vitro, usually fewer than ten divisions. These cells retain their diploid karyotype, i.e., they have the chromosome number and morphology of their tissues of origin. They also retain some of the differentiated characteristics that they possessed in vivo. Because of this, these cells support the replication of a wide range of viruses. Primary cultures derived from monkey kidneys, mouse fetuses, and chick embryos are commonly used for diagnostic purposes and laboratory experiments.

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Some primary cells can be passed through secondary and several subsequent subcultures while retaining their original morphological characteristics and karyotype. Subcultures will have fewer cell types than primary cultures. After 20 to 50 passages in vitro, these diploid cell strains usually undergo a crisis in which their growth rate slows and they eventually die out. Diploid strains of

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Inverted phase contrast microscopes are used for visualizing the cells. Before using the microscope, check that the phase rings are aligned. Inverted phase contrast microscopes are used for visualizing the cells. Before using the microscope, check that the phase rings are aligned. الإستعمال.

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## 9.3 Preservation and storage of tissue cells / الحفظ و التخزين

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## MAINTENANCE

Cultures should be examined daily, observing the morphology, the color of the medium and the density of the cells. A tissue culture log should be maintained that is separate from your regular laboratory notebook. The log should contain: the name of the cell line, the medium components and any alterations to the standard medium, the dates on which the cells were split and/or fed, a calculation of the doubling time of the culture (this should be done at least once during the semester), etc.

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### 9.4 Harvesting and refeeding culture cells / الحصاد

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Trypsin, collagenase, or pronase, usually in combination with EDTA, causes cells to detach from the growth surface. This method is fast and reliable but can damage the cell surface by digesting exposed cell surface proteins. The proteolysis reaction can be quickly terminated by the addition of complete medium containing serum.

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EDTA alone can also be used to detach cells and seems to be gentler on the cells than trypsin.

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The standard procedure for detaching adherent cells is as follows:

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- humidity is required
- gas phase - bicarbonate conc. and CO<sub>2</sub> tension in equilibrium
- visible light - can have an adverse effect on cells; cells should be cultured in the dark and exposed to room light as little as possible

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The viability of cells can be observed visually using an inverted phase contrast microscope. Live cells are phase bright; suspension cells are typically rounded and somewhat symmetrical; adherent cells will form projections when they attach to the

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growth surface. Viability can also be assessed using the vital dye, trypan blue, which is excluded by live cells but accumulates in dead cells. Cell numbers are determined using a Hemacytometer.

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Assume all cultures are hazardous since they may harbor latent viruses or other organisms that are uncharacterized.

The following safety precautions should also be observed:

- pipetting: use pipette aids to prevent ingestion and keep aerosols down to a minimum
- no eating, drinking, or smoking
- wash hands after handling cultures and before leaving the lab
- decontaminate work surfaces with disinfectant (before and after)
- autoclave all waste
- use safety cabinet when working with hazardous organisms.
- use aseptic technique
- dispose of all liquid waste after each experiment.

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The cells should be monitored daily for morphology and growth characteristics, fed every 2 to 3 days, and subcultured when necessary. A minimum of two 25 cm<sup>2</sup> flasks should be carried for each cell line. Each time the cells are subcultured, a viable cell count should be done, the subculture dilutions should be noted, and, after several passages, a doubling time determined. As soon as you have enough cells, several vials should be frozen away and stored in liquid N<sub>2</sub>. One vial from each freeze down

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should be thawed 1-2 weeks after freezing to check for viability. These frozen stocks will prove to be vital if any of your cultures become contaminated.

the minimal nutritional requirements of cultured cells :mixture of salts, amino acids, vitamins and cofactors, carbohydrates, and horse serum. By eliminating one component at a time , then determined which nutrients were essential for cell growth. His minimum essential medium (MEM) contains 13 amino acids (human tissue in vivo requires only eight), 8 vitamins and cofactors, glucose as an energy source, and a physiological salt solution that is isotonic to the cell. The pH is maintained at 7.2 to 7.4 by  $\text{NaHCO}_3$  in equilibrium with  $\text{CO}_2$ . The pH indicator phenol red is usually incorporated into the medium; it turns red-purple if the medium is basic, yellow if the medium is acidic, and remains red-orange if the pH is in the right range. Serum in concentrations of 1 to 10% must be added to the medium to provide the cells with additional poorly defined factors, without which most cells will not grow. Most mammalian cells are incubated at  $37^\circ\text{C}$ ; avian, reptilian, and arthropod cells may grow best at higher or lower temperatures.

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When adherent cells become semi-confluent, subculture using 2 mM EDTA or trypsin/EDTA.

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- Remove medium from culture dish and wash cells in a balanced salt solution without  $\text{Ca}^{++}$  or  $\text{Mg}^{++}$ . Remove the wash solution.
- Add enough trypsin-EDTA

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ice-cold freezing medium: 90% calf serum/10% DMSO or glycerol with ethylene glycol , at  $10^6 - 10^7$  cells/ml. Keep cells on ice.

- Transfer 1 ml aliquots to freezer vials on ice.
- Place in a Mr. Frosty container that is at room temperature and that has sufficient isopropanol.
- Place the Mr. Frosty in the -70°C freezer overnight. Note: Cells should be exposed to freezing medium for as little time as possible prior to freezing
- Next day, transfer to liquid nitrogen (DON'T FORGET).

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## 9.7.6 Viable cell counts with Hemacytometer / حساب الخلايا

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USE A HEMACYTOMETER TO DETERMINE TOTAL CELL COUNTS AND VIABLE CELL NUMBERS.

Blue is one of several stains recommended for use in dye exclusion procedures for viable cell counting. This method is based on the principle that live cells do not take up certain dyes, whereas dead cells do.

- Prepare a cell suspension, either directly from a cell culture or from a concentrated or diluted suspension and combine 20  $\mu\text{l}$  of cells with 20  $\mu\text{l}$  of trypan blue suspension (0.4%). Mix thoroughly and allow to stand for 5-15 minutes.
- With the cover slip in place, transfer a small amount of trypan blue-cell suspension to both chambers of the Hemacytometer by carefully touching the edge of the cover slip with the pipette tip and allowing each chamber to fill by capillary action. Do not overfill or underfill

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Add distilled H<sub>2</sub>O to 1 litre. Filter sterilize.

Saline A is usually prepared as a 10X solution and stored at -20 °C

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
Add distilled H<sub>2</sub>O to 1 liter. Filter sterilize. Minimum essential medium Eagle is usually purchased as a preweighed mixture or as a sterile solution. The above formulation has Earle's salts.

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- Incubator
- Centrifuge
- Forceps and scissors
- Alcohol

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### 10.3 Protocol

1-Desinfect the surface of the egg over the air sac with scissors or forceps, break the shell.

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2-Sterilize forceps by dipping in alcohol and flaming. Cool forceps, then peel away the shell over the air sac

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3-Sterilize forceps again, and pull back the shell membrane and chorioallantoic membrane to expose the embryo.

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4-Resterilize the forceps, grasp the embryo loosely around the neck, and remove the entire embryo from the egg to a sterile Petri dish.

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5-Using forceps + scissors decapitate and eviscerate the embryo. Mince the embryo carcass into very small fragment with scissors .

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6- Add about 10ml sterile saline A to tissue fragment in the Petri dish,

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Swirl gently for 1 to 2 minutes, and carefully pour the entire content into a 125 ml Erlenmeyer flask ,tilt flask ,for decant saline A. Discard the saline A.

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7-Add 10 ml of sterile warm trypsin and pour into centrifuge filtrate at 1000 rpm for 10 minutes ,discard supernatant and resuspend the cells in 20 ml growth medium .

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8-Mix well for counting in a hemacytometer.

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N.B: be sure to keep your cell sterile.

9-Dilute 0.1ml of cell suspension with 0.9ml of MEM and count cell in a hemacytometer

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10-Adjust concentration to  $5 \times 10^5$  cells per ml growth medium and dispense into cultures tubes.

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11-Incubate at  $37^\circ\text{C}$  until monolayered (6-7 days).

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**Note:** be sure examine cell culture each day ,if the medium is yellow may need to be adjust by  $7.5\% \text{NaHCO}_3$  .if floating cell are present or the color of the medium indicate a basic PH the medium should be changed.

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